

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21

**LUMI-CELL[®] ER ASSAY
AGONIST PROTOCOL**

**National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative
Toxicological Methods (NICEATM)**

**Developed by:
Xenobiotic Detection Systems, Inc.
1601 E. Geer St., Suite S
Durham, NC 27704
*12 March 2009***

22 **TABLE OF CONTENTS**

23 **LIST OF ACRONYMS AND ABBREVIATIONS vi**

24 **LIST OF FIGURES vii**

25 **LIST OF TABLES viii**

26 **1.0 Purpose.....1**

27 **2.0 Sponsor1**

28 **2.1 Substance Inventory and Distribution Management3**

29 **3.0 Definitions.....3**

30 **4.0 Testing Facility and Key Personnel.....4**

31 4.1 Testing Facility4

32 4.2 Key Personnel4

33 **5.0 Identification of Test and Control Substances4**

34 5.1 Test Substances4

35 5.2 Controls.....4

36 **6.0 Overview of General Procedures For Agonist Testing.....5**

37 6.1 Range Finder Testing.....6

38 6.2 Comprehensive Testing6

39 **7.0 Materials for LUMI-CELL® ER Agonist Testing7**

40 7.1 BG1Luc4E2 Cells7

41 7.2 Technical Equipment7

42 7.3 Reference Standard, Controls, and Tissue Culture Supplies9

43 **8.0 Preparation of Tissue Culture Media and Solutions10**

44 8.1 RPMI 1640 Growth Medium (RPMI)10

45 8.2 Estrogen-free DMEM Medium.....11

46	8.3	1X Trypsin Solution.....	11
47	8.4	1X Lysis Solution	12
48	8.5	Reconstituted Luciferase Reagent.....	12
49	9.0	Overview of Propagation and Experimental Plating of	
50		BG1Luc4E2 Cells	13
51	9.1	Procedures for Thawing Cells and Establishing Tissue Cultures	13
52	9.1.1	Thawing Cells	13
53	9.1.2	Establishing Tissue Cultures.....	14
54	9.2	Ongoing Tissue Culture Maintenance, Conditioning in Estrogen-free	
55		Medium, and Plating Cells for Experimentation	16
56	9.2.1	Ongoing Tissue Culture Maintenance	17
57	9.2.2	Conditioning in Estrogen-free Medium	17
58	9.2.3	Plating Cells Grown in Estrogen-free DMEM for	
59		Experimentation.....	18
60	10.0	Preparation of Test Substances	21
61	10.1	Determination of Test Substance Solubility	21
62	10.2	Preparation of Reference Standards, Control and Test Substances	22
63	10.2.1	Preparation of Reference Standard and Positive Control Stock	
64		Solutions	22
65	10.2.1.1	E2 Stock Solution	22
66	10.2.1.2	Methoxychlor Stock Solution	22
67	10.2.2	Preparation of Reference Standard and DMSO Control for	
68		Range Finder Testing.....	23
69	10.2.2.1	Preparation of E2 Reference Standard for	
70		Range Finder Testing.....	23

71	10.2.2.2	Preparation of DMSO Control for Range	
72		Finder Testing	24
73	10.2.3	Preparation of Test Substance Dosing Solutions for Range	
74		Finder Testing	24
75	10.2.4	Preparation of Reference Standard and Positive Control	
76		Dosing Solutions for Comprehensive Testing	25
77	10.2.4.1	Preparation of E2 Reference Standard for	
78		Comprehensive Testing	25
79	10.2.4.2	Preparation of Methoxychlor Control Dosing	
80		Solution for Comprehensive Testing	26
81	10.2.4.3	Preparation of DMSO Control Dosing	
82		Solution for Comprehensive Testing	26
83	10.2.5	Preparation of Test Substance Dosing Solutions for	
84		Comprehensive Testing	26
85	10.2.5.1	Preparation of Test Substance 1:2 Serial	
86		Dilutions for Comprehensive Testing	27
87	10.2.5.2	Preparation of Test Substance 1:5 Serial	
88		Dilutions for Comprehensive Testing	28
89	11.0	General Procedures for the Testing of Coded Substances	28
90	11.1	Application of Reference Standard, Controls, and Test Substances	29
91	11.1.1	Preparation of Excel® Data Analysis Template for Range Finder	
92		Testing	29
93	11.1.2	Preparation of Excel® Data Analysis Template for Comprehensive	
94		Testing	30
95	11.2	Visual Evaluation of Cell Viability	31
96	11.3	Lysis of Cells for LUMI-CELL® ER	32
97	11.4	Measurement of Luminescence	32

98	11.5	Data Analysis	32
99	11.5.1	Collection and Adjustment of Luminometer Data for	
100		Range Finder Testing	33
101	11.5.2	Collection and Adjustment of Luminometer Data for	
102		Comprehensive Testing	35
103	11.5.3	Determination of Outliers	37
104	11.5.4	Acceptance Criteria.....	38
105	11.5.4.1	Range Finder Testing.....	38
106	11.5.4.2	Comprehensive Testing	38
107	12.0	Range Finder Testing	39
108	13.0	Comprehensive Testing	44
109	14.0	Compilation of the Historical Quality Control Database	46
110	14.1	Compilation of Historical Quality Control Database for	
111		Phase IIa Testing.....	46
112	14.1.1	DMSO Control.....	46
113	15.0	Quality Testing of Materials	46
114	15.1	Tissue Culture Media	47
115	15.2	G418.....	48
116	15.3	DMSO.....	48
117	15.4	Plastic Tissue Culture Materials	49
118	16.0	References.....	50
119			

119	LIST OF ACRONYMS AND ABBREVIATIONS	
120	13 mm test tube	13 x 100 mm glass test tubes
121	DMEM	Dulbecco's Modification of Eagle's Medium
122	DMSO	Dimethyl Sulfoxide
123	DMSO control	1% v/v dilution of DMSO in tissue culture media used as a
124		vehicle control
125	E2	17 β -estradiol
126	E2 reference standard	11 Point Serial Dilution of 17 β -estradiol reference standard
127		for the LUMI-CELL® ER agonist assay
128	EC ₅₀ value	Concentration that produces a half-maximal response as
129		calculated using the four parameter Hill function.
130	ER	Estrogen Receptor
131	Estrogen-free DMEM	DMEM (phenol red free) supplemented with 1%
132		Penicillin/Streptomycin, 2% L-Glutamine, and 5%
133		Charcoal-dextran treated FBS
134	FBS	Fetal Bovine Serum
135	G418	Gentamycin
136	Methoxychlor	<i>p,p'</i> -Methoxychlor
137	Methoxychlor control	3.13 μ g/mL Methoxychlor Weak Positive Control for the
138		LUMI-CELL® ER Agonist Assay
139	RPMI	RPMI 1640 growth medium
140	TA	Transcriptional Activation
141	T25	25 cm ² tissue culture flask
142	T75	75 cm ² tissue culture flask
143	T150	150 cm ² tissue culture flask
144		

144

LIST OF FIGURES

145 Figure 7-1 pGudLuc7.ERE Plasmid7

146 Figure 9-1 Hemocytometer Counting Grid.....19

147 Figure 12-1 Agonist Range Finder Test Plate Layout39

148 Figure 12-2 Agonist Range Finder (example 1).....41

149 Figure 12-3 Agonist Range Finder (example 2).....42

150 Figure 12-4 Agonist Range Finder (example 3).....42

151 Figure 12-5 Agonist Range Finder (example 4).....43

152 Figure 12-6 Agonist Range Finder (example 5).....43

153 Figure 12-7 Agonist Range Finder (example 6).....44

154 Figure 13-1 Agonist Comprehensive Test Plate Layout.....45

155

156

156 **LIST OF TABLES**

157 Table 6-1 Concentrations of E2 Reference Standard Used in
158 Comprehensive Testing5

159 Table 10-1 Preparation of E2 Stock Solution.....22

160 Table 10-2 Preparation of E2 Reference Standard Dosing Solution for Range
161 Finder Testing24

162 Table 10-3 Preparation of Test Substance Dosing Solutions for Range
163 Finder Testing24

164 Table 10-4 Preparation of E2 Reference Standard Dosing Solution for
165 Comprehensive Testing26

166 Table 10-5 Preparation of Test Substance 1:2 Serial Dilutions for
167 Comprehensive Testing27

168 Table 10-6 Preparation of Test Substance 1:5 Serial Dilutions for
169 Comprehensive Testing28

170 Table 11-1 Visual Observation Scoring32

171 Table 11-2 Q Test Values.....37

172

173

174 **1.0 PURPOSE**

175 This protocol is designed to evaluate coded test substances for potential estrogen receptor (ER)
176 agonist activity using the LUMI-CELL® ER assay.

177 **2.0 SPONSOR**

178 The National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative
179 Toxicological Methods (NICEATM), P.O. Box 12233 Research Triangle Park, NC 27709

180

181 William S. Stokes, DVM, DACLAM
182 Rear Admiral, U.S. Public Health Service
183 Chief Veterinary Officer, USPHS
184 Director, NICEATM
185 National Institute of Environmental Health Sciences, NIH, DHHS
186 Bldg. 4401, Room 3129, MD EC-14
187 79 T.W. Alexander Drive
188 Research Triangle Park, NC 27709
189 Phone: 919-541-7997
190 Fax: 919-541-0947
191 Email: stokes@niehs.nih.gov

192

193 Raymond Tice, Ph.D.
194 Deputy Director, NICEATM
195 National Institute of Environmental Health Sciences
196 MD EC-17, P.O. Box 12233
197 Research Triangle Park, NC 27709
198 Phone: 919-541-4482
199 FAX: 919-541-0947
200 Email: tice@niehs.nih.gov

201

202

202 David Allen, Ph.D.
203 Principal Investigator
204 ILS, Inc./Contractor supporting NICEATM
205 National Institute of Environmental Health Sciences
206 MD EC-17, P.O. Box 12233
207 Research Triangle Park, NC 27709
208 Phone: 919-316-4664
209 FAX: 919-541-0947
210 Email: allen7@niehs.nih.gov
211
212 Frank Deal, M.S.
213 Staff Toxicologist
214 ILS, Inc./Contractor supporting NICEATM
215 National Institute of Environmental Health Sciences
216 MD EC-17, P.O. Box 12233
217 Research Triangle Park, NC 27709
218 Phone: 919-316-4587
219 FAX: 919-541-0947
220 Email: dealf@niehs.nih.gov
221
222 Patricia Ceger, M.S.
223 Project Coordinator/Technical Writer
224 ILS, Inc./Contractor supporting NICEATM
225 National Institute of Environmental Health Sciences
226 MD EC-17, P.O. Box 12233
227 Research Triangle Park, NC 27709
228 Phone: 919-316-4556
229 Fax: 919-541-0947
230 E-Mail: cegerp@niehs.nih.gov
231
232

232 **2.1 Substance Inventory and Distribution Management**

233 Cynthia Smith, Ph.D.

234 Chemistry Resources Group Leader

235 National Institute of Environmental Health Sciences

236 MD EC-06, P.O. Box 12233

237 Research Triangle Park, NC 27709

238 Phone: 919-541-3473

239

240 **3.0 DEFINITIONS**

241 • **Dosing Solution:** The test substance, control substance, or reference standard
242 solution, which is to be placed into the tissue culture wells for experimentation.

243 • **Raw Data:** Raw data includes information that has been collected but not
244 formatted or analyzed, and consists of the following:

245 ○ Data recorded in the Study Notebook

246 ○ Computer printout of initial luminometer data

247 ○ Other data collected as part of GLP compliance, e.g.:

248 ■ Equipment logs and calibration records

249 ■ Test substance and tissue culture media preparation logs

250 ■ Cryogenic freezer inventory logs

251 • **Soluble:** Test substance exists in a clear solution without visible cloudiness or
252 precipitate.

253 • **Study Notebook:** The study notebook contains recordings of all activities related
254 to the conduct of the LUMI-CELL® ER agonist assay.

255 • **Test Substances:** Substances supplied to the testing laboratories that are coded
256 and distributed such that only the Project Officer, Study Management Team
257 (SMT), and the Substance Inventory and Distribution Management have
258 knowledge of their true identity. The test substances will be purchased, aliquoted,

259 coded, and distributed by the Supplier under the guidance of the NIEHS/NTP
260 Project Officer and the SMT.

261 **4.0 TESTING FACILITY AND KEY PERSONNEL¹**

262 **4.1 Testing Facility**

263 Xenobiotic Detection Systems, Inc. (XDS), 1601 E. Geer St., Suite S, Durham, NC 27704

264 **4.2 Key Personnel**

- 265 • Study Director: John Gordon, Ph.D.
- 266 • Quality Assurance Director: Mr. Andrew

267 **5.0 IDENTIFICATION OF TEST AND CONTROL SUBSTANCES**

268 **5.1 Test Substances**

269 Test substances are coded and will be provided to participating laboratories by the Substance
270 Inventory and Distribution Management team.

271 **5.2 Controls**

272 Controls for the ER agonist protocol are as follows:

273 *Vehicle control (dimethyl sulfoxide [DMSO]):* 1% (v/v) DMSO (CASRN 67-68-5) diluted in
274 tissue culture media.

275 *Reference standard (17 β -estradiol [E2]):* Three concentrations of E2 (CASRN 50-28-2) in
276 duplicate for range finder testing and a serial dilution consisting of 11 concentrations of E2 in
277 duplicate for comprehensive testing

278 *Positive control (p,p'-Methoxychlor [methoxychlor]):* Methoxychlor (CASRN 72-43-5), 3.13
279 μ g/mL in tissue culture media, used as a weak positive control.

¹ Testing facility and personnel information are provided as an example.

280 **6.0 OVERVIEW OF GENERAL PROCEDURES FOR AGONIST TESTING**

281 All experimental procedures are to be carried out under aseptic conditions and all solutions,
 282 glassware, plastic ware, pipettes, etc., shall be sterile. All methods and procedures shall be
 283 documented in the study notebook.

284 Agonist range finder testing is conducted on 96-well plates using four concentrations of E2
 285 (5.00×10^{-5} , 1.25×10^{-5} , 3.13×10^{-6} and 7.83×10^{-7} $\mu\text{g/mL}$) in duplicate as the reference standard
 286 and four replicate wells for the DMSO control. Range finder testing uses all wells of the 96-well
 287 plate to test six substances as seven point 1:10 serial dilutions in duplicate.

288 Comprehensive testing is conducted on 96-well plates using 11 concentrations of E2 in duplicate
 289 as the reference standard (**Table 6-1**). Four replicate wells for the DMSO control and four
 290 replicate wells for the methoxychlor control are included on each plate. Comprehensive testing
 291 uses all wells of the 96-well plate to test 2 substances as 11 point serial dilutions in triplicate.

292 **Table 6-1 Concentrations of E2 Reference Standard Used in Comprehensive Testing**

E2 Concentrations ¹		
1.00×10^{-4}	6.25×10^{-6}	3.92×10^{-7}
5.00×10^{-5}	3.13×10^{-6}	1.95×10^{-7}
2.50×10^{-5}	1.56×10^{-6}	9.78×10^{-8}
1.25×10^{-5}	7.83×10^{-7}	

293 ¹Concentrations are presented in $\mu\text{g/mL}$.

294 Visual observations for cell viability are conducted for all experimental plates just prior to
 295 luminescence measurements, as outlined in **Section 11.2**.

296 Luminescence data, measured in relative light units (RLUs), is corrected for background
 297 luminescence by subtracting the mean RLU value of the vehicle control (DMSO) wells from the
 298 RLU measurements for each of the other wells of the 96-well plate. Data is then transferred into
 299 Excel® data management spreadsheets and GraphPad PRISM® 4.0 statistical software, graphed,
 300 and evaluated as follows:

- 301 • A response is considered positive for agonist activity when the average adjusted
 302 RLU for a given concentration is greater than the mean RLU value plus three
 303 times the standard deviation for the vehicle control.

304 • Any response below this threshold is considered negative for agonist activity.

305 For substances that are positive at one or more concentrations, the concentration that causes a
306 half-maximal response (EC₅₀) is calculated using a Hill function analysis. The Hill function is a
307 four-parameter logistic mathematical model relating the substance concentration to the response
308 (typically following a sigmoidal curve) using the equation below:

$$309 \quad Y = \text{Bottom} + \frac{\text{Top} - \text{Bottom}}{1 + 10^{(\log \text{EC}_{50} - X) \text{HillSlope}}}$$

310 where Y = response (i.e., relative light units); X = the logarithm of concentration; Bottom = the
311 minimum response; Top = the maximum response; log EC₅₀ = the logarithm of X as the response
312 midway between Top and Bottom; and HillSlope describes the steepness of the curve. The model
313 calculates the best fit for the Top, Bottom, HillSlope, and EC₅₀ parameters. See **Section 11.6.5**
314 for more details.

315 Acceptance or rejection of a test is based on evaluation of reference standard and control results
316 from each experiment conducted on a 96-well plate. Results for these controls are compared to
317 historical results compiled in the historical database, as seen in **Section 14.0**.

318 **6.1 Range Finder Testing**

319 Agonist range finding for coded substances consists of a seven point, 1:10 serial dilution using
320 duplicate wells per concentration. Concentrations for comprehensive testing are selected based
321 on the response observed in range finder testing. If necessary, a second range finder test can be
322 conducted to clarify the optimal concentration range to test (see **Section 12.0**).

323 **6.2 Comprehensive Testing**

324 Comprehensive agonist testing for coded substances consists of 11 point, serial dilutions, with
325 each concentration tested in triplicate wells of the 96-well plate. Three separate experiments are
326 conducted for comprehensive testing on three separate days, except during Phases III and IV of
327 the validation effort, in which comprehensive testing experiments are conducted once (see
328 **Section 13.0**).


329 7.0 MATERIALS FOR LUMI-CELL® ER AGONIST TESTING

330 This section provides the materials needed to conduct LUMI-CELL® ER testing, with associated
331 brand names/vendors² in brackets.

332 7.1 BG1Luc4E2 Cells:

333 Human ovarian cancer cell line stably transfected with a plasmid containing an estrogen response
334 element pGudLuc7.0 (**Figure 7-1**) [XDS].

335 **Figure 7-1 pGudLuc7.ERE Plasmid.**



336

337 7.2 Technical Equipment:

338 All technical equipment may be obtained from Fisher Scientific International, Inc. (Liberty Lane
339 Hampton, NH, USA 03842). Equivalent technical equipment from another commercial source
340 can be used.

- 341
- Analytical balance (Cat. No. 01-910-320)
 - 342
 - Berthold Orion 1 Microplate Luminometer [Berthold CatNo.: Orion 1 MPL3] or
343 equivalent and dedicated computer
 - 344
 - Biological safety hood, class II, and stand (Cat. No. 16-108-99)

²Brand names and vendors should not be considered an endorsement by the U.S. Government or any member of the U.S. Government; such information is provided as examples.

- 345 • Centrifuge (low speed, tabletop with swinging bucket rotor) (Cat. No. 04-978-50
- 346 centrifuge, and 05-103B rotor)
- 347 • Combustion test kit (CO₂ monitoring) (Cat. No. 10-884-1)
- 348 • Drummond diaphragm pipetter (Cat. No. 13-681-15)
- 349 • Freezers, –20°C (Cat. No. 13-986-150), and –70°C (Cat. No. 13-990-86)
- 350 • Hand tally counter (Cat. No. 07905-6)
- 351 • Hemocytometer, cell counter (Cat. No. 02-671-5)
- 352 • Light microscope, inverted (Cat. No. 12-561-INV)
- 353 • Light microscope, upright (Cat. No. 12-561-3M)
- 354 • Liquid nitrogen flask (Cat. No. 11-675-92)
- 355 • Micropipetter, repeating (Cat. No. 21-380-9)
- 356 • Pipettors, air displacement, single channel (0.5 –10µl (Cat. No. 21-377-191), 2 –
- 357 20 µl (Cat. No. 21-377-287), 20 – 200 µl (Cat. No. 21-377-298), 200 - 1000 µl
- 358 (Cat. No. 21-377-195))
- 359 • Refrigerator/freezer (Cat. No. 13-986-106A)
- 360 • Shaker for 96-well plates (Cat. No. 14-271-9)
- 361 • Sodium hydroxide (Cat. No. 5318-500)
- 362 • Sonicating water bath (Cat. No. 15-335-30)
- 363 • Tissue culture incubator with CO₂ and temperature control (Cat. No. 11-689-4)
- 364 • Vacuum pump with liquid trap (side arm Erlenmeyer) (Cat. No. 01-092-29)
- 365 • Vortex mixer (Cat. No. 12-814)

366 Equipment should be maintained and calibrated as per GLP guidelines and individual laboratory
367 SOPs.

368

368 7.3 Reference Standard, Controls, and Tissue Culture Supplies

369 All tissue culture reagents must be labeled to indicate source, identity, storage conditions and
370 expiration dates. Tissue culture solutions must be labeled to indicate concentration, stability
371 (where known), and preparation and expiration dates.

372 Equivalent tissue culture media and sera from another commercial source can be used, but must
373 first be tested as described in **Section 15.0** to determine suitability for use in this test method.

374 The following are the necessary tissue culture reagents and possible commercial sources (in
375 brackets) based on their use in the pre-validation studies:

- 376 • BackSeal-96/384, white adhesive bottom seal for 96-well and 384-well microplate
377 [Perkin-Elmer, Cat. No. 6005199]
- 378 • 17 β -estradiol (CAS RN: 50-28-2) [Sigma-Aldrich, Cat. No. E8875]
- 379 • Cryovial, 2 mL (Corning Costar) [Fisher Scientific Cat. No. 03-374-21]
- 380 • Culture tube 13 x 100mm (case) [Thomas Scientific Cat. No.: 10009186R38]³
- 381 • Culture tube, 50 mL conical (Corning Costar) [Fisher Scientific Cat. No. 05-
382 526C]
- 383 • DMSO, U.S.P. analytical grade. [Sigma-Aldrich, Cat. No. 34869-100ML]
- 384 • Dulbecco's Modification of Eagle's Medium (DMEM), containing 4.5 g/L
385 glucose, with sodium pyruvate, without phenol red or L-glutamine
386 [Mediatech/Cellgro, Cat. No. 17-205-CV]
- 387 • Fetal Bovine Serum [Mediatech/Cellgro Cat. No. MT 35-010-CV]
- 388 • Fetal Bovine Serum, charcoal/dextran treated, triple 0.1 μ m sterile filtered
389 [Hyclone, Cat. No. SH30068.03]
- 390 • Gentamycin Sulfate (G418), 50 mg/mL [Mediatech/Cellgro Cat. No. 30-234-CR]
- 391 • L-glutamine, 29.2 mg/mL [Cellgro, Cat. No. 25005-CI]

³ If glass tubes can not be obtained from Thomas Scientific, the preference is for flint glass, then lime glass, then borosilicate glass.

- 392 • Luciferase Assay System (10-Pack) [Promega Cat. No. E1501]
- 393 • Lysis Solution 5X [Promega, Cat. No. E1531]
- 394 • Methoxychlor (CAS RN: 72-43-5) [Sigma-Aldrich, Cat. No. 49054]
- 395 • Penicillin/streptomycin solution, 5000 I.U. penicillin, 5000 µg/mL streptomycin
- 396 [Cellgro, Cat. No. 30-001-CI].
- 397 • Phosphate buffered saline (PBS, 1X) without calcium and magnesium [Cellgro,
- 398 Cat. No. 21-040-CV]
- 399 • Pipettes, serological: 2.0 mL [Sigma-Aldrich, Cat. No. P1736], 5.0 mL [Sigma-
- 400 Aldrich, Cat. No. P1986], 25 mL [Sigma-Aldrich, Cat. No. P2486]
- 401 • RPMI 1640 medium, containing L-glutamine [Mediatech, Cat. No. 10-040-CV]
- 402 • Tissue culture flasks (Corning-Costar): 25 cm² (T25) [Fisher Cat. No. 10-126-28];
- 403 75 cm² (T75) [Fisher Cat. No. 10-126-37]; and 150 cm² (T150) [Fisher Cat. No.
- 404 10-126-34]
- 405 • Tissue culture plates (Corning-Costar): 96-well [Thomas Scientific Cat. No.
- 406 6916A05]
- 407 • Trypsin (10X), 2.5% in Hank's balanced salt solution (HBSS), without calcium
- 408 and magnesium, without phenol red [Cellgro, Cat. No. 25-054-CI].

409 All reagent lot numbers and expiration dates must be recorded in the study notebook.

410 **8.0 PREPARATION OF TISSUE CULTURE MEDIA AND SOLUTIONS**

411 All tissue culture media and media supplements must be quality tested before use in experiment
412 (see Section 15.0).

413 **8.1 RPMI 1640 Growth Medium (RPMI)**

414 RPMI 1640 is supplemented with 0.9% Pen-Strep and 8.0% FBS to make RPMI growth medium
415 (RPMI).

416 Procedure for one 549 mL bottle:

- 417 1. Remove FBS from -70°C freezer, and Pen-Strep from -20°C freezer and allow to
418 equilibrate to room temperature.
- 419 2. Add 44 mL of FBS and 5 mL Pen-Strep to the bottle of RPMI 1640.
- 420 3. Label RPMI bottle as indicated in **Section 7.3**

421 *Store at 2-8°C for no longer than six months or until the shortest expiration date of any media*
422 *component.*

423 **8.2 Estrogen-Free DMEM Medium**

424 DMEM is supplemented to contain 4.5% charcoal/dextran treated FBS, 1.9% L-glutamine, 0.9%
425 Pen-Strep.

426 Procedure for one 539 mL bottle:

- 427 1. Remove charcoal/dextran treated FBS from -70°C freezer, and L-glutamine and
428 Pen-Strep from -20°C freezer and allow to equilibrate to room temperature.
- 429 2. Add 24 mL of charcoal/dextran treated FBS, 10 mL L-glutamine, and 5 mL Pen-
430 Strep to one 500 mL bottle of DMEM.
- 431 3. Label estrogen-free DMEM bottle as indicated in **Section 7.3**

432 *Store at 2-8°C for no longer than six months or until the shortest expiration date of any media*
433 *component..*

434 **8.3 1X Trypsin Solution**

435 1X Trypsin solution is prepared by dilution from a 10X premixed stock solution. The 10X stock
436 solution should be stored in 10 mL aliquots in a -20°C freezer.

437 Procedure for making 100 mL of 1X trypsin:

- 438 1. Remove a 10 mL aliquot of 10X trypsin from -20°C freezer and allow to
439 equilibrate to room temperature.
- 440 2. Aliquot 1 mL Trypsin (10X) along with 9 mL of 1X PBS into ten 15 mL sterile
441 centrifuge tubes.
- 442 3. Label 1X trypsin aliquots as indicated in **Section 7.3**

443 *1X Trypsin should be stored at -20°C.*

444 **8.4 1X Lysis Solution**

445 Lysis solution is prepared by dilution from a 5X premixed stock solution. Both the 5X and 1X
446 solutions can be repeatedly freeze-thawed.

447 The procedure for making 10 mL of 1X lysis solution:

- 448 1. Thaw the 5X Promega Lysis solution and allow it to reach room temperature.
- 449 2. Remove 2 mL of 5X solution and place it in a 15 mL conical centrifuge tube.
- 450 3. Add 8 mL of distilled, de-ionized water to the conical tube.
- 451 4. Cap and shake gently until solutions are mixed.

452 *Store at -20°C for no longer than 1 year from receipt.*

453 **8.5 Reconstituted Luciferase Reagent**

454 Luciferase reagent consists of two components, luciferase buffer and lyophilized luciferase
455 substrate.

456 For long term storage, unopened containers of the luciferase buffer and lyophilized luciferase
457 substrate can be stored at -70°C for up to one year.

458 To reconstitute luciferase reagent:

- 459 1. Remove luciferase buffer and luciferase substrate from -70°C freezer and allow
460 them to equilibrate to room temperature.
- 461 2. Add 10 mL of luciferase buffer solution to luciferase substrate container and swirl
462 or vortex gently to mix; the Luciferase substrate should readily go into solution.
- 463 3. After solutions are mixed, aliquot to a 15mL centrifuge tube.
- 464 4. Store complete solution at -20°C.

465 Reconstituted luciferase reagent is stable for up to 1 month at -20°C.

466 **9.0 OVERVIEW OF PROPAGATION AND EXPERIMENTAL PLATING OF** 467 **BG1Luc4E2 CELLS**

468 The BG1Luc4E2 (BG-1) cells are stored in liquid nitrogen in 2 mL cryovials. BG-1 cells are
469 grown as a monolayer in tissue culture flasks in a dedicated tissue culture incubator at 37°C ±
470 1°C, 90% ± 5% humidity, and 5.0% ± 1% CO₂/air. The cells should be examined, on a daily
471 basis during working days, under an inverted phase contrast microscope and any changes in
472 morphology and/or adhesive properties must be noted in the study notebook.

473 Two T150 flasks containing cells at 80 to 90% confluence will usually yield a sufficient number
474 of cells to fill three 96-well plates for use in experiments.

475 **9.1 Procedures for Thawing Cells and Establishing Tissue Cultures**

476 Warm all of the tissue culture media and solutions to room temperature by placing them under
477 the tissue culture hood several hours before use.

478 All tissue culture media, media supplements, and tissue culture plasticware must be quality
479 tested before use in experiments (**Section 15.0**).

480 9.1.1 Thawing Cells

- 481 1. Remove a cryovial of frozen BG-1 cells from the liquid nitrogen flask.
- 482 2. Facilitate rapid thawing by loosening the top slightly (do not remove top) to
483 release trapped gasses and retightening it. Roll vial between palms.
- 484 3. Use a micropipette to transfer cells to a 50 mL conical centrifuge tube.
- 485 4. Rinse cryovial twice with 1X PBS and add PBS rinse material to the conical tube.
- 486 5. Add 20 mL of RPMI to the conical tube.
- 487 6. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
488 for an additional 5 minutes.
- 489 7. Aspirate media from pellet and re-suspend it in 5 mL RPMI, drawing the pellet
490 repeatedly through a 1.0 mL serological pipette to break up any clumps of cells.
- 491 8. Transfer cells to a T25 flask, place them in an incubator (see conditions in
492 **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).

493 9.1.2 Establishing Tissue Cultures

494 Once cells have reached 80 to 90% confluence, transfer the cells to a T75 flask by performing,
495 for example, the following steps:

- 496 1. Remove the T25 flask from the incubator.
- 497 2. Aspirate the RPMI, then add 5 mL 1X PBS, making sure that the cells are coated
498 with PBS.
- 499 3. Aspirate 1X PBS, then add 1 to 2 mL 1X trypsin to the T25 flask, gently swirling
500 the flask to coat all cells with the trypsin.
- 501 4. Place the flask in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 502 5. Detach cells by hitting the side of the flask sharply against the palm or heel of the
503 hand.
- 504 6. Confirm cell detachment by examination under an inverted microscope. If cells
505 have not detached, return the flask to the incubator for an additional 2 minutes,
506 then hit the flask again.
- 507 7. After cells have detached, add 5 mL PBS, and transfer the suspended cells to a 50
508 mL centrifuge tube. Wash the flask one additional time with 5 mL PBS.
- 509 8. Immediately add 20 mL RPMI to the conical tube to inhibit further cellular
510 digestion by residual trypsin.
- 511 9. Pellet the cells by centrifugation, as described in **Section 9.1.1**, and re-suspend the
512 cells in 10 mL RPMI medium.
- 513 10. Draw the pellet repeatedly through a 25 mL serological pipette to break up
514 clumps of cells
- 515 11. Transfer cells to a T75 flask, then place the flask in an incubator (see conditions
516 in **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).

517 When cells have reached 80% to 90% confluency, transfer them into a T150 flask by performing,
518 for example, the following steps:

- 519 12. Remove the T75 flask from the incubator, aspirate the old media and add 5 mL
520 1X PBS.
- 521 13. Aspirate 1X PBS, add 2 mL of 1X trypsin to the flask, and place it in an incubator
522 (see conditions in **Section 9.0**) for 5 to 10 min.
- 523 14. Repeat steps **5** through **11** in **Section 9.1.2**, re-suspending the pellet in 20 mL of
524 RPMI.
- 525 15. Transfer cells to a T150 flask and place it in the incubator (see conditions in
526 **Section 9.0**) and grow to 80% to 90% confluence (approximately 48 to 72 hrs).
- 527 16. Remove the T150 flask from the incubator.
- 528 17. Aspirate the RPMI and add 5 mL 1X PBS.
- 529 18. Aspirate 1X PBS and add 3 mL 1X trypsin to the T150 flask, making sure that the
530 cells are coated with the trypsin.
- 531 19. Incubate cells in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 532 20. Detach cells by hitting the side of the flask sharply against the palm or heel of the
533 hand.
- 534 21. Confirm cell detachment by examination under an inverted microscope. If cells
535 have not detached, return the flask to the incubator for an additional 2 minutes,
536 then hit the flask again.
- 537 22. After cells have detached, add 5 mL 1X PBS and transfer the suspended cells
538 from the T150 flask to a 50 mL conical tube. Add an additional 5 mL PBS to the
539 flask, swirl around the flask, and then transfer the PBS to the 50 mL conical tube.
- 540 23. Immediately add 20 mL RPMI to the conical tube to inhibit further cellular
541 digestion by residual trypsin.
- 542 24. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
543 for an additional 5 minutes.

- 544 25. Aspirate the media from the pellet and re-suspend it in 40 mL RPMI, drawing the
545 pellet repeatedly through a 25 mL serological pipette to break up any clumps of
546 cells.
- 547 26. Transfer 20 mL of cell suspension to each of two T150 flasks, place them in an
548 incubator (see conditions in **Section 9.0**) and grow to 80% to 90% confluence
549 (approximately 48 to 72 hrs).

550 **9.2 Ongoing Tissue Culture Maintenance, Conditioning in Estrogen-free Medium,**
551 **and Plating Cells for Experimentation**

552 The following procedure is used to condition the BG1Luc4E2 cells to an estrogen-free
553 environment prior to plating the cells in 96-well plates for analysis of estrogen dependent
554 induction of luciferase activity.

555 To start the tissue culture maintenance and estrogen-free conditioning, split the two T150 culture
556 flasks into four T150 flasks. Two of these flasks will be used for continuing tissue culture and
557 will use the RPMI media mentioned above. The other two flasks will be cultured in estrogen-free
558 DMEM for experimental use. Extra care must be taken to avoid contaminating the estrogen-free
559 cells with RPMI.

- 560 1. Remove both T150 flasks from the incubator.
- 561 2. Aspirate the medium and rinse the cells with 5 mL 1X PBS.
- 562 3. Aspirate 1X PBS, then add 3 mL 1X trypsin to the flasks, gently swirling the flask
563 to coat all cells with the trypsin.
- 564 4. Incubate cells in the incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 565 5. Detach cells by hitting the side of the flask sharply against the palm or heel of the
566 hand.
- 567 6. Confirm cell detachment by examination under an inverted microscope. If cells
568 have not detached, return the flask to the incubator for an additional 2 minutes,
569 then hit the flask again.
- 570 7. After cells have detached, add 5 mL 1X PBS to the first T150 flask and transfer
571 the suspended cells to the second T150 flask.

- 572 8. Transfer the contents of both flasks to a 50 mL conical tube. Repeat step 7 with an
573 additional 5 mL 1X PBS and transfer to the 50 mL conical tube.
- 574 9. Immediately add 20 mL estrogen-free DMEM to the 50 mL conical tube to inhibit
575 further cellular digestion by residual trypsin.
- 576 10. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
577 for an additional 5 minutes.
- 578 11. Aspirate media from pellet and re-suspend it in 4 mL estrogen-free DMEM,
579 drawing the pellet repeatedly through a 1 mL serological pipette to break up
580 clumps of cells.

581 At this point, cells are ready to be divided into the ongoing tissue culture and estrogen-free
582 conditioning groups.

583 9.2.1 Ongoing Tissue Culture Maintenance

- 584 1. Add 20 mL RPMI to two T150 flasks.
- 585 2. Add 220 µl G418 to the RPMI in the T150 flasks
- 586 3. Add 1 mL of cell suspension from **9.2 step 11** to each flask.
- 587 4. Place T150 flasks in tissue culture incubator (see conditions in **Section 9.0**) and
588 grow to 80% to 90% confluence (approximately 48 to 72 hrs).
- 589 5. Tissue culture medium may need to be changed 24 hours after addition of G418 to
590 remove cells that have died because they do not express reporter plasmid.
- 591 6. G418 does not need to be added to the flasks a second time.
- 592 7. Repeat **Section 9.2 steps 1-11** for ongoing tissue culture maintenance.

593 9.2.2 Conditioning in Estrogen-free Medium

- 594 1. Add 20 mL estrogen-free DMEM to two T150 flasks.
- 595 2. Add 150 µL G418 to the estrogen-free DMEM in the T150 flasks.
- 596 3. Add 1 mL of cell suspension from **Section 9.2 step 11** to each flask.


- 597 4. Tissue culture medium may need to be changed 24 hours after addition of G418 to
598 remove cells that have died because they do not express reporter plasmid.
- 599 5. G418 does not need to be added to the flasks a second time.
- 600 6. Place the T150 flasks in the incubator (see conditions in **Section 9.0**) and grow to
601 80% to 90% confluence (approximately 48 to 72 hrs).

602 9.2.3 Plating Cells Grown in Estrogen-free DMEM for Experimentation

- 603 1. Remove the T150 flasks that have been conditioned in estrogen-free DMEM for
604 48 to 72 hours from the incubator.
- 605 2. Aspirate the medium, then rinse the cells with 5 mL 1X PBS.
- 606 3. Aspirate 1X PBS, then add 3 mL 1X trypsin to the flasks, gently swirling the flask
607 to coat all cells with the trypsin.
- 608 4. Place the flasks in an incubator (see conditions in **Section 9.0**) for 5 to 10 min.
- 609 5. Detach cells by hitting the side of the flask sharply against the palm or the heel of
610 the hand.
- 611 6. Confirm cell detachment by examination under an inverted microscope. If cells
612 have not detached, return the flask to the incubator for 2 additional minutes, then
613 hit the flask again.
- 614 7. After cells have detached, add 5 mL 1X PBS and transfer the suspended cells
615 from the T150 flask to a 50 mL conical tube. Add an additional 5 mL PBS to the
616 flask, gently swirl around the flask, and then transfer to the 50 mL conical tube.
- 617 8. Immediately add 20 mL estrogen-free DMEM to each conical tube to inhibit
618 further cellular digestion by residual trypsin.
- 619 9. Centrifuge at 1000 x g for eight min. If a pellet of cells has not formed, centrifuge
620 for an additional 5 minutes.
- 621 10. Aspirate the media from the pellet and re-suspend it in 20 mL DMEM, drawing
622 the pellet repeatedly through a 25 mL serological pipette to break up any clumps
623 of cells.

- 624 11. Pipette 15 µL of the cell suspension into the “v” shaped slot on the
 625 hemocytometer. Ensure that the solution covers the entire surface area of the
 626 hemocytometer grid, and allow cells to settle before counting.
- 627 12. Using 100x magnification, view the counting grid.
- 628 13. The counting grid on the hemocytometer consists of nine sections, four of which
 629 are counted (upper left, upper right, lower left, and lower right, see **Figure 9-1**).
 630 Each section counted consists of four by four grids. Starting at the top left and
 631 moving clockwise, count all cells in each of the four by four grids. Some cells
 632 will be touching the outside borders of the square, but only count those that touch
 633 the top and right borders of the square. This value is then used in the calculation
 634 below to get to the desired concentration of 200,000 cells/mL.
- 635

Figure 9-1 Hemocytometer Counting Grid.



636 The volume of each square is 10^{-4} mL, therefore:

637
$$\text{Cells/mL} = (\text{average number per grid}) \times 10^{-4} \text{ mL} \times 1/(\text{starting dilution}).$$

638 Starting dilution: 20 mL (for T150 flasks)

639

640 Harvested cells for a T150 flask are suspended in 20 mL of estrogen-free DMEM and sampled
 641 for determination of concentration of cells/mL.

642

643 Example Calculation:

644 • Grids 1, 2, 3, and 4 are counted and provide the following data:

645 ○ 50, 51, 49, and 50: average number of cells per grid is equal to 50.

646 Cells/mL = 50 cells per grid ÷ 10^{-4} mL volume of grid = 50×10^4 cells/mL (or 500,000
647 cells/mL)

648 Total # of Cells Harvested = 500,000 cells/mL x 20 mL

649 Desired Concentration (or Concentration_{Final}) = 200,000 cells/mL

650 Formula: (Concentration_{Final} x Volume_{Final} = Concentration_{Initial} x Volume_{Initial})

651 Concentration_{Final} = 200,000 cells/mL

652 Concentration_{Initial} = 500,000 cells/mL

653 Volume_{Initial} = 20 mL

654 Volume_{Final} – to be solved for.

655 Therefore: 200,000 cells/mL x Volume_{Final} = 500,000 cells/mL x 20 mL

656 Solving for Volume_{Final} we find = 50 mL

657 Therefore, add 30 mL of estrogen-free DMEM to the cell suspension for a total volume of 50
658 mL, which will yield the desired concentration of 200,000 cells/mL for plating.

659 14. This dilution scheme will give a concentration of 200,000 cells/mL. 200 µL of
660 this cell suspension is used for each well of a 96-well plate (i.e., 40,000 cells per
661 well).

662 15. Remove a 96-well plate from its sterile packaging. Use a repeater pipetter to
663 pipette 200 µL of cell suspension into each well for to be used for the testing of
664 coded substances, reference standard and controls (**note:** add 200 µL of estrogen-
665 free DMEM only to any wells not being used for testing).

666 16. Incubate plate(s) in an incubator (see conditions in **Section 9.0**) for a minimum of
667 24 hours, but no longer than 48 hours before dosing.

668 Two T150 flasks containing cells at 80% to 90% confluence will typically yield sufficient cells
669 to fill four 96-well plates.

670 **10.0 PREPARATION OF TEST SUBSTANCES**

671 The solvent used for dissolution of test substances is 100% DMSO. All test substances should be
672 allowed to equilibrate to room temperature before being dissolved and diluted. Test substance
673 solutions (except for reference standards and controls) should not be prepared in bulk for use in
674 subsequent tests. Test substances are to be used within 24 hours of preparation. Solutions should
675 not have noticeable precipitate or cloudiness.

676 All information on weighing, solubility testing, and calculation of final concentrations for test
677 substances, reference standards and controls is to be recorded in the study notebook.

678 **10.1 Determination of Test Substance Solubility**

- 679 1. Prepare a 100 mg/mL solution of the test substance in 100% DMSO in a 4 mL
680 conical tube.
- 681 2. Vortex to mix.
- 682 3. If the test substance does not dissolve at 100 mg/mL, prepare a 10 mg/mL
683 solution and vortex as above.
- 684 4. If the test substance does not dissolve at 10 mg/mL solution, prepare a 1 mg/mL
685 solution in a 4 mL conical tube and vortex as above.
- 686 5. If the test substance does not dissolve at 1 mg/mL, prepare a 0.1 mg/mL solution
687 in a 4 mL conical tube and vortex as above.
- 688 6. Continue testing, using 1/10 less substance in each subsequent attempt until test
689 substance is solubilized in DMSO.

690 Once the test substance has fully dissolved in 100% DMSO, the test substance is ready to be
691 used for LUMI-CELL® ER testing.

692 The Testing Facility shall forward the results from the solubility tests assay to the SMT through
693 the designated contacts in electronic format and hard copy upon completion of testing.

694 **10.2 Preparation of Reference Standards, Control and Test Substances**

695 All “dosing solutions” of test substance concentrations are to be expressed as µg/mL in the study
696 notebook and in all laboratory reports.

697 All information on preparation of test substances, reference standards and controls is to be
698 recorded in the study notebook.

699 **10.2.1 Preparation of Reference Standard and Positive Control Stock Solutions**

700 Stock solutions of E2 and methoxychlor are prepared in 100% DMSO and stored at room
701 temperature for up to three years or until the expiration date listed in the certificate of analysis
702 for that substance.

703 **10.2.1.1 E2 Stock Solution**

704 The final concentration of the E2 stock solution is 1.0×10^{-2} µg/mL. Prepare the E2 stock as
705 shown in **Table 10-1**.

706 **Table 10-1 Preparation of E2 Stock Solution**

Step #	Action	DMSO	E2 Concentration
1	Make a 10 mg/mL stock solution in 100% DMSO in a 4mL vial.	-	10 mg/mL
2	Transfer 10 µL E2 solution from Step #1 to a new 4 mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	100 µg/mL
3	Transfer 10 µL E2 solution from Step #2 to a new 4mL vial.	Add 990 µL of 100% DMSO. Vortex to mix.	1 µg/mL
4	Transfer 10 µL E2 solution from Step #3 to a 13 mm test tube to create the working solution.	Add 990 µL of 100% DMSO. Vortex to mix.	1.0×10^{-2} µg/mL

707

708 **10.2.1.2 Methoxychlor Stock Solution**

709 The final concentration of the methoxychlor stock solution is 313 µg/mL.

710 To prepare the methoxychlor stock solution, proceed as follows:

711 1. Make a 10 mg/mL stock solution of Methoxychlor in 100% DMSO in a 4 mL
712 vial.

713 2. Remove 94 µL of the methoxychlor solution and place it in a new 4 mL vial.

714 3. Add 2.906 mL of 100% DMSO to the 4mL vial and gently vortex to mix.

715 10.2.2 Preparation of Reference Standard and DMSO Control for Range Finder Testing

716 Range finder testing is conducted on 96-well plates using four concentrations of E2 in duplicate
717 as the reference standard. Four replicate wells are used for the DMSO control. All wells on the
718 96 well plate are used during range finder testing.

719 Store dosing solutions at room temperature. Use within 24 hours of preparation.

720 10.2.2.1 *Preparation of E2 Reference Standard for Range Finder Testing*

721 To make E2 dosing solutions:

722 1. label four 4 mL conical tubes with numbers 1 through 4 and place them in a tube
723 rack

724 2. label four 13 mm glass test tubes with numbers 1 through 4, place them in a tube
725 rack and add 600 µL of estrogen-free DMEM to each tube

726 Prepare dilutions to give final concentrations of the E2 as shown in **Table 10-2**.

727

727 **Table 10-2 Preparation of E2 Reference Standard Dosing Solution for Range Finder**
 728 **Testing**

Tube Number	100% DMSO	E2 ¹	Estrogen-free DMEM ²	Final Volume	E2 Concentration
1	6 µL	6 µl of 1.0 x 10 ⁻² µg/mL working solution	600 µL	606 µL	5.00 x 10 ⁻⁵ µL
2	18 µL	6 µL of 1.0 x 10 ⁻² µg/mL working solution	600 µL	606 µL	1.25 x 10 ⁻⁵ µL
3	18 µL	6 µL from conical tube #2	600 µL	606 µL	3.13 x 10 ⁻⁶ µL
4	18 µL	6 µL from conical tube #3	600 µL	606 µL	7.83 x 10 ⁻⁷ µL

729 ¹Add specified volume of 100% DMSO and 6 µl of the specified E2 solution to labeled 4 mL conical
 730 tubes, and vortex.

731 ²Transfer 6 µL of DMSO/E2 solution from 4 mL conical tube to labeled 13 mL glass tubes containing
 732 DMEM and vortex.

733 10.2.2.2 *Preparation of DMSO Control for Range Finder Testing*

734 1. Add 10 µL of 100% DMSO to a 13 mm glass tube.

735 2. Add 1000 µL of estrogen-free DMEM to tube and vortex.

736 10.2.3 Preparation of Test Substance Dosing Solutions for Range Finder Testing

737 Range finder experiments are used to determine the concentrations of test substance to be used
 738 during comprehensive testing. Agonist range finding for coded substances consists of seven
 739 point, 1:10 serial dilutions run in duplicate.

740 To make dosing solutions for coded substances:

741 1. label seven 4 mL conical tubes with numbers 1 through 7 and place them in a tube
 742 rack

743 2. label seven 13 mm glass test tubes with numbers 1 through 7, place them in a tube
 744 rack and add 600 µL of estrogen-free DMEM to each tube

745 Prepare dilutions as shown in **Table 10-3**.

746 **Table 10-3 Preparation of Test Substance Dosing Solutions for Range Finder Testing**

Tube Number	100% DMSO	Test Substance ¹	Transfer	Estrogen-free DMEM	Final Volume
1	-	6 µL of test substance solution from Section 10.1 step 10	6 µL	600 µL	606 µL

Tube Number	100% DMSO	Test Substance ¹	Transfer	Estrogen-free DMEM	Final Volume
2	90 µL	10 µL of test substance solution from Section 10.1 step 10	6 µL	600 µL	606 µL
3	90 µL	10 µL from conical tube #2	6 µL	600 µL	606 µL
4	90 µL	10 µL from conical tube #3	6 µL	600 µL	606 µL
5	90 µL	10 µL from conical tube #4	6 µL	600 µL	606 µL
6	90 µL	10 µL from conical tube #5	6 µL	600 µL	606 µL
7	90 µL	10 µL from conical tube #6	6 µL	600 µL	606 µL

747 ¹Add specified volume of 100% DMSO and test substance solution to labeled 4 mL conical tubes, and
 748 vortex.

749 ²Transfer 6 µL of DMSO/E2 solution from 4 mL conical tube to labeled 13 mL glass tubes containing
 750 DMEM and vortex.

751
 752 Determination of whether a substance is positive in range finder testing and selection of starting
 753 concentrations for comprehensive testing will be discussed in **Section 12.0**.

754 10.2.4 Preparation of Reference Standard and Positive Control Dosing Solutions for 755 Comprehensive Testing

756 Comprehensive testing is conducted on 96-well plates using 11 concentrations of E2 in duplicate
 757 as the reference standard. Four replicate wells for the DMSO control and three replicate wells for
 758 the methoxychlor control are included on each plate.

759 Store dosing solutions at room temperature. Use within 24 hours of preparation.

760 10.2.4.1 *Preparation of E2 Reference Standard for Comprehensive Testing*

761 To make E2 dosing solutions:

- 762 1. label 11 4 mL conical tubes with numbers 1 through 11 and place them in a tube
 763 rack
- 764 2. label 11 13 mm glass test tubes with numbers 1 through 11, place them in a tube
 765 rack and add 600 µL of DMEM to each tube

766 Prepare dilutions to give final concentrations of E2 as shown in **Table 10-4**.

767

767 **Table 10-4 Preparation of E2 Reference Standard Dosing Solution for**
 768 **Comprehensive Testing**

Tube Number	100% DMSO	E2 ¹	Estrogen-free DMEM ²	Final Volume	E2 Concentration
1	-	6 µl of 1.0 x 10 ⁻² µg/mL working solution	600 µL	606 µL	1.00 x 10 ⁻⁴ µL
2	6 µL	6 µL of 1.0 x 10 ⁻² µg/mL working solution	600 µL	606 µL	5.00 x 10 ⁻⁵ µL
3	6 µL	6 µL from conical tube #2	600 µL	606 µL	2.50 x 10 ⁻⁵ µL
4	6 µL	6 µL from conical tube #3	600 µL	606 µL	1.25 x 10 ⁻⁵ µL
5	6 µL	6 µL from conical tube #4	600 µL	606 µL	6.25 x 10 ⁻⁶ µL
6	6 µL	6 µL from conical tube #5	600 µL	606 µL	3.13 x 10 ⁻⁶ µL
7	6 µL	6 µL from conical tube #6	600 µL	606 µL	1.56 x 10 ⁻⁶ µL
8	6 µL	6 µL from conical tube #7	600 µL	606 µL	7.83 x 10 ⁻⁷ µL
9	6 µL	6 µL from conical tube #8	600 µL	606 µL	3.92 x 10 ⁻⁷ µL
10	6 µL	6 µL from conical tube #9	600 µL	606 µL	1.95 x 10 ⁻⁷ µL
11	6 µL	6 µL from conical tube #10	600 µL	606 µL	9.78 x 10 ⁻⁸ µL

769 ¹Add specified volume of 100% DMSO and 6 µl of the specified E2 solution to labeled 4 mL conical
 770 tubes, and vortex.

771 ²Transfer 6 µL of DMSO/E2 solution from 4 mL conical tube to labeled 13 mL glass tubes containing
 772 DMEM and vortex.

773

774 10.2.4.2 *Preparation of Methoxychlor Control Dosing Solution for Comprehensive Testing*

775 1. Add 10 µL of the 313 µg/mL methoxychlor to a 13 mm glass tube.

776 2. Add 1000 µL of estrogen-free DMEM to the tube and vortex.

777 10.2.4.3 *Preparation of DMSO Control Dosing Solution for Comprehensive Testing*

778 1. Add 10 µL of 100% DMSO to four 13 mm tubes (solvent/negative controls).

779 2. Add 1000 µL of estrogen-free DMEM to the tube and vortex.

780 10.2.5 Preparation of Test Substance Dosing Solutions for Comprehensive Testing

781 Comprehensive testing experiments are used to determine whether a substance possesses ER

782 agonist activity in the LUMI-CELL® ER test method. Agonist comprehensive testing for coded

783 substances consists of either an 11 point 1:2 serial dilution or an 11 point 1:5 serial dilution,

784 depending on the results from range finder testing (see **Section 12.0**) with each concentration
 785 tested in triplicate wells of the 96-well plate.

786 10.2.5.1 *Preparation of Test Substance 1:2 Serial Dilutions for*
 787 *Comprehensive Testing*

788 Start the 11-point serial dilution according to criteria in **Section 12.0**.

789 To make test substance 1:2 serial dilutions for comprehensive testing:

- 790 1. label eleven 4 mL conical tubes with numbers 1 through 11 and place them in a
 791 tube rack
- 792 2. label eleven 13 mm glass test tubes with numbers 1 through 11, place them in a
 793 tube rack and add 800 µL of estrogen-free DMEM to each tube

794 Prepare dilution of test substance as shown in **Table 10-6**.

795 **Table 10-5 Preparation of Test Substance 1:2 Serial Dilutions for Comprehensive**
 796 **Testing**

Tube Number	100% DMSO	Test Substance ¹	Transfer	Estrogen-free DMEM	Final Volume
1	-	8 µL of highest concentration of test substance solution	8 µL	800 µL	808 µL
2	8 µL	8 µL of highest concentration of test substance solution	8 µL	800 µL	808 µL
3	8 µL	8 µL from conical tube #2	8 µL	800 µL	808 µL
4	8 µL	8 µL from conical tube #3	8 µL	800 µL	808 µL
5	8 µL	8 µL from conical tube #4	8 µL	800 µL	808 µL
6	8 µL	8 µL from conical tube #5	8 µL	800 µL	808 µL
7	8 µL	8 µL from conical tube #6	8 µL	800 µL	808 µL
8	8 µL	8 µL from conical tube #7	8 µL	800 µL	808 µL
9	8 µL	8 µL from conical tube #8	8 µL	800 µL	808 µL
10	8 µL	8 µL from conical tube #9	8 µL	800 µL	808 µL
11	8 µL	8 µL from conical tube #10	8 µL	800 µL	808 µL

797 ¹Add specified volume of 100% DMSO and test substance solution to labeled 4 mL conical tubes, and
 798 vortex.

799

800 10.2.5.2 *Preparation of Test Substance 1:5 Serial Dilutions for Comprehensive*
 801 *Testing*

802 Start the 11-point serial dilution according to criteria in **Section 12.0**.

803 To make test substance 1:5 serial dilutions for comprehensive testing:

804 3. label eleven 4 mL conical tubes with numbers 1 through 11 and place them in a
 805 tube rack

806 4. label eleven 13 mm glass test tubes with numbers 1 through 11, place them in a
 807 tube rack and add 800 µL of estrogen-free DMEM to each tube

808 Prepare dilution of test substance as shown in **Table 10-6**.

809 **Table 10-6 Preparation of Test Substance 1:5 Serial Dilutions for Comprehensive**
 810 **Testing**

Tube Number	100% DMSO	Test Substance ¹	Transfer	Estrogen-free DMEM	Final Volume
1	-	8 µL of highest concentration of test substance solution	8 µL	800 µL	808 µL
2	16 µL	4 µL of highest concentration of test substance solution	8 µL	800 µL	808 µL
3	16 µL	4 µL from conical tube #2	8 µL	800 µL	808 µL
4	16 µL	4 µL from conical tube #3	8 µL	800 µL	808 µL
5	16 µL	4 µL from conical tube #4	8 µL	800 µL	808 µL
6	16 µL	4 µL from conical tube #5	8 µL	800 µL	808 µL
7	16 µL	4 µL from conical tube #6	8 µL	800 µL	808 µL
8	16 µL	4 µL from conical tube #7	8 µL	800 µL	808 µL
9	16 µL	4 µL from conical tube #8	8 µL	800 µL	808 µL
10	16 µL	4 µL from conical tube #9	8 µL	800 µL	808 µL
11	16 µL	4 µL from conical tube #10	8 µL	800 µL	808 µL

811 ¹Add specified volume of 100% DMSO and test substance solution to labeled 4 mL conical tubes, and
 812 vortex.
 813

814 **11.0 GENERAL PROCEDURES FOR THE TESTING OF CODED SUBSTANCES**

815 Range finder experiments are used to determine the concentrations of test substance to be used
 816 during comprehensive testing. Comprehensive testing experiments are used to determine whether
 817 a substance possesses ER agonist activity in the LUMI-CELL® ER assay.

818 General procedures for range finder and comprehensive are similar. For specific details (such as
819 plate layout) of range finder testing see **Section 12.0**. For specific details of comprehensive
820 testing, see **Section 13.0**.

821 **11.1 Application of Reference Standard, Controls, and Test Substances**

- 822 1. Remove seeded 96-well plates from the incubator, inspect them using an inverted
823 microscope. Only use plates in which the cells in all wells giving a score of 1
824 according to **Table 11-1**.
- 825 2. Remove medium by inverting the plate onto blotter paper. Gently tap plate against
826 the bench surface to remove residual liquid trapped in the wells.
- 827 3. Add 200 µL of reference standard, control, or test substance to each well (see
828 **Sections 12.0** and **13.0** for specific plate layouts).
- 829 4. Return plates to incubator and incubate (see **Section 9.0** for details) for 19 to 24
830 hours to allow maximal induction of luciferase activity in the cells.

831 **11.1.1 Preparation of Excel® Data Analysis Template For Range Finder Testing**

- 832 1. In Excel®, open a new “AgRFTemplate” and save it with the appropriate project
833 name as indicated in the NICEATM Style Guide.
- 834 2. Fill out the table at the top of the “Raw Data” worksheet with information
835 regarding the Microplate reader used, Reading Direction, No. of Intervals, Tot.
836 Meas. Time/Well (s), etc. (**note**: this information can be permanently added to the
837 default template “AgRFTemplate” on a laboratory specific basis).
- 838 3. Add the following information regarding the assay to the “Compound Tracking”
839 worksheet.
 - 840 ▪ Plate # - Enter the experiment ID or plate number into cell E1
 - 841 ▪ Cell Lot # - Enter the passage or lot number of the cells used for this
842 experiment into cell B5
 - 843 ▪ DMSO and Media Lot #'s – Enter the lot numbers for the DMSO and
844 Media in cells B6 and B7

- 845 ▪ Test Substance Code – Enter the test substance codes into cells C13 to
846 C18
- 847 ▪ Name: Enter the experimenter name into cell G6
- 848 ▪ Date: Enter the experiment date in the format day\month\year into cell
849 G10
- 850 ▪ Comments: - Enter any comments about the experiment in this box (e.g.,
851 plate contaminated)
- 852 4. Enter the following substance testing information to the “List” page:
- 853 ▪ Concentration – Type in the test substance concentration in µg/ml in
854 descending order.
- 855 ▪ Also add any replicate-specific comments on this page (e.g, spilled tube,
856 etc.), in the comments section
- 857 ▪ All of the remaining cells on the List tab should populate automatically.
- 858 ▪ The “Template”, “Compound Mixing” and “Visual Inspection” tabs
859 should automatically populate with the information entered into the
860 Compound Tracking and List tabs.
- 861 5. Save the newly named project file.
- 862 6. Print out either the “List” or “Template” page for help with dosing the 96-well
863 plate. Sign and date the print out and store in study notebook.

864 11.1.2 Preparation of Excel® Data Analysis Template for Comprehensive Testing

- 865 1. In Excel®, open a new “AgCTTemplate” and save it with the appropriate project
866 name as indicated in the NICEATM Style Guide.
- 867 2. Fill out the table at the top of the “Raw Data” worksheet with information
868 regarding the Microplate reader used, Reading Direction, No. of Intervals, Tot.
869 Meas. Time/Well (s), etc. (**note**: this information can be permanently added to the
870 default template “AgCTTemplate” on a laboratory specific basis).

871

- 872 3. On the “Compound Tracking” tab, enter the following information:
- 873 ▪ Plate # - Enter the experiment ID or plate number into cell E1
- 874 ▪ Cell Lot # - Enter the passage or lot number of the cells used for this
- 875 experiment into cell C5
- 876 ▪ DMSO and Media Lot #'s – Enter the lot numbers for the DMSO and
- 877 Media in cells C6 and C7
- 878 ▪ Test Substance Code – Enter the test substance codes into cells C15 and
- 879 C16. Enter the test substance dilution into cells E25 and E26.
- 880 ▪ Name: Enter the experimenter name into cell G6
- 881 ▪ Date: Enter the experiment date in the format day\month\year into cell
- 882 G10
- 883 ▪ Comments: - Enter any comments about the experiment in this box (e.g.,
- 884 plate contaminated)
- 885 4. Enter substance testing concentrations to the “List” page. Also add any replicate-
- 886 specific comments on this page (e.g, spilled tube, etc.).
- 887 5. Save the newly named project file.
- 888 6. Print out either the “List” or “Template” page for help with dosing the 96-well
- 889 plate. Sign and date the print out and store in study notebook.

890 11.2 Visual Evaluation of Cell Viability

- 891 1. 19 to 24 hours after dosing the plate, remove the plate from the incubator and
- 892 remove the media from the wells by inverting the plate onto blotter paper. Gently
- 893 tap plate against the bench surface to remove residual liquid trapped in the wells.
- 894 2. Use a repeat pipetter to add 50 µL 1X PBS to all wells. Immediately remove PBS
- 895 by inversion.
- 896 3. Using an inverted microscope, inspect all of the wells used in the 96-well plate
- 897 and record the visual observations using the scores in **Table 11-1**.

898 **Table 11-1 Visual Observation Scoring**

Viability Score	Brief Description ¹
1	Normal Cell Morphology and Cell Density
2	Altered Cell Morphology and/or Small Gaps between Cells
3	Altered Cell Morphology and/or Large Gaps between Cells
4	Few (or no) Visible Cells
P	Wells containing precipitation are to be noted with “P”

899 ¹Reference photomicrographs are provided in the LUMI-CELL® ER Validation Study “Visual Observation Cell
900 Viability Manual.”
901

902 **11.3 Lysis of Cells for LUMI-CELL® ER**

- 903 1. Apply the reflective white backing tape to the bottom of the 96-well plate (this
904 will increase the effectiveness of the luminometer).
- 905 2. Add 30 µL 1X lysis reagent to the assay wells and place the 96-well plate on an
906 orbital shaker for one minute.
- 907 3. Remove plate from shaker and measure luminescence (as described in **Section**
908 **11.4**).

909 **11.4 Measurement of Luminescence**

910 Luminescence is measured in the range of 300 to 650 nm, using an injecting luminometer and
911 with software that controls the injection volume and measurement interval. Light emission from
912 each well is expressed as RLU per well. The luminometer output is saved as raw data in an
913 Excel® spreadsheet. A hard copy of the luminometer raw data should be signed, dated and stored
914 in the study notebook.

915 **11.5 Data Analysis**

916 LUMI-CELL® ER uses an Excel® spreadsheet to collect and adjust the RLU values obtained
917 from the luminometer and a GraphPad Prism® template to analyze and graph data. The Excel®
918 spreadsheet subtracts background luminescence (average DMSO solvent control RLU value)
919 from test substance, reference standard and control RLU values. Plate induction is calculated
920 using these corrected RLU values. Test substance, reference standard, and control RLU values
921 are then adjusted relative to the highest E2 reference standard RLU value, which is set to 10,000.
922 After adjustment, values are transferred to GraphPad Prism® for data analysis and graphing.

923 11.5.1 Collection and Adjustment of Luminometer Data for Range Finder Testing

924 The following steps describe the procedures required to populate the Excel® spreadsheet that has
925 been configured to collect and adjust the RLU values obtained from the luminometer.

- 926 1. Open the raw data file and the corresponding experimental Excel® spreadsheet
927 from **Section 11.1.1**.
- 928 2. Copy the raw data using the Excel® copy function, then paste the copied data into
929 cell B19 of the “RAW DATA” tab in the experimental Excel® spreadsheet using
930 the **Paste Special – Values** command. This position corresponds to position A1 in
931 the table labeled Table 1 in this tab.
- 932 3. Examine the DMSO data in Table 1 of the Excel® spreadsheet to determine
933 whether there are any potential outliers. See **Section 11.6.2** for further explanation
934 of outlier determinations.
- 935 4. If an outlier is identified, perform the following steps to remove the outlier from
936 calculations:
- 937 ▪ correct the equation used to calculate DMSO background in Table 1 [e.g.,
938 if outlier is located in cell F26, adjust the calculation in cell H40 to read
939 =AVERAGE(G26:I26)]
 - 940 ▪ then correct the equation used to calculate the average DMSO value in
941 Table 2 [e.g., following the above example, adjust cell M42 to read
942 =AVERAGE(G26:I26)]
 - 943 ▪ then correct the equation used to calculate the standard deviation of the
944 DMSO value in Table 2 [e.g., following the above example, adjust cell
945 M43 to read =STDEV(G36:I36)]
- 946 5. Excel® will automatically subtract the background (the average DMSO control
947 value) from all of the RLU values in Table 1 and populate Table 2 with these
948 adjusted values.
- 949 6. To calculate plate induction, identify the cell containing the E2a replicate in Table
950 1, plate row H that has the highest RLU value (i.e., cell B26, C26, D26, or E26).

- 951 7. Click into cell D14 and enter the cell number from the previous step into the
952 numerator.
- 953 8. Identify the cell containing the E2b replicate in Table 1, plate row H that has the
954 highest RLU value (i.e., cell J26, K26, L26, or M26).
- 955 9. Click into cell E14 and enter the cell number from the previous step into the
956 numerator.
- 957 10. Click on the “ER Agonist Report” worksheet.
- 958 11. The data for the E2 reference standard, methoxychlor, and DMSO replicates
959 populate the left portion (columns A – F) of the spreadsheet. The data is
960 automatically placed in an Excel® graph.
- 961 12. To set the highest RLU value for the reference standard to 10,000 RLU, go to cell
962 D2 of “ER Agonist Report” tab and check the formula contained within that cell.
963 The divisor should be the cell number of the cell containing the highest Mean E2
964 RLU value ((i.e., cell A16, A17, A18, or A19).
- 965 13. Open the “Visual Observation Scoring” worksheet. Enter the visual observation
966 scores for each well on the 96-well plate. This data will be linked to the “ER
967 Agonist Report” worksheet.
- 968 14. After the testing results have been evaluated and reviewed for quality control,
969 enter the following information into the Compound Tracking worksheet:
- 970 ▪ Enter pass/fail results for plate reference standard and control parameters
971 into the Plate Pass/Fail Table
 - 972 ▪ Enter information from the testing of coded substances into the Testing
973 Results Table
 - 974 ▪ Reviewer Name – Enter the name of the person who Reviewed\QC’ed the
975 data into cell A34
 - 976 ▪ Date – Enter the date on which the data was reviewed into cell D34
977

977 11.5.2 Collection and Adjustment of Luminometer Data for Comprehensive Testing

978 The following steps describe the procedures required to populate the Excel® spreadsheet that has
979 been configured to collect and adjust the RLU values obtained from the luminometer.

- 980 1. Open the raw data file and the corresponding experimental Excel® spreadsheet
981 from **Section 11.1.2**.
- 982 2. Copy the raw data using the Excel® copy function, then paste the copied data into
983 cell B16 of the “RAW DATA” worksheet in the experimental Excel® spreadsheet
984 using the **Paste Special – Values** command. This position corresponds to position
985 A1 in the table labeled Table 1 in this worksheet.
- 986 3. Fill out the table at the top of the “Raw Data” worksheet with information
987 regarding the Microplate reader used, Reading Direction, No. of Intervals, Tot.
988 Meas. Time/Well (s), etc. If desired, this information can be added to the
989 Laboratory Template File.
- 990 4. Examine the DMSO data in Table 1 of the Excel® spreadsheet to determine
991 whether there are any potential outliers. See **Section 11.6.2** for further explanation
992 of outlier determinations.
- 993 5. If an outlier is identified, perform the following steps to remove the outlier from
994 calculations:
- 995 ▪ correct the equation used to calculate DMSO background in Table 1 [e.g.,
996 if outlier is located in cell M17, adjust the calculation in cell H37 to read
997 =AVERAGE(M16,M18:M19)]
 - 998 ▪ then correct the equation used to calculate the DMSO mean and SD
999 values [e.g., following the above example, adjust cell M39 to read
1000 =AVERAGE(M28,M30:M31), and adjust cell M40 to read
1001 =STDEV(M28,M30:M31)]
- 1002 6. Excel® will automatically subtract the background (the average DMSO control
1003 value) from all of the RLU values in Table 1 and populate Table 2 with these
1004 adjusted values.

- 1005 7. To calculate plate induction, identify the cell in containing the E2 replicate in
1006 Table 1, plate row G that has the highest RLU value.
- 1007 8. Click into cell D11 and enter the cell number from the previous step into the
1008 numerator.
- 1009 9. Identify the cell containing the E2 replicate in plate row H that has the highest
1010 RLU value.
- 1011 10. Click into cell E11 and enter the cell number from the previous step into the
1012 numerator.
- 1013 11. Open the “ER Agonist Report” worksheet.
- 1014 12. The data for the E2 reference standard, methoxychlor, and DMSO replicates
1015 populate the left portion (columns A – E) of the spreadsheet. The data is
1016 automatically placed in an Excel® graph.
- 1017 13. To set the highest RLU value for the reference standard to 10,000 RLU, go to cell
1018 E2 of “ER Agonist Report” tab and check the formula contained within that cell.
1019 The divisor should be the cell number of the cell containing the highest Avg E2
1020 RLU value (cells A16 through A26).
- 1021 14. Open the “Visual Observation Scoring” worksheet. Enter the visual observation
1022 scores for each well on the 96-well plate. This data will be linked to the “ER
1023 Agonist Report” worksheet.
- 1024 15. Copy the data from the “ER Agonist Report” worksheet into GraphPad Prism® for
1025 the calculation of EC₅₀ values and to graph experimental results as indicated in the
1026 NICEATM Prism® Users Guide.
- 1027 16. After the testing results have been evaluated and reviewed for quality control,
1028 enter the following information into the Compound Tracking worksheet:
- 1029 ▪ Enter pass/fail results for plate reference standard and control parameters
1030 into the Plate Pass/Fail Table
 - 1031 ▪ Enter information from the testing of coded substances into the Testing
1032 Results Table

- 1033 ▪ Reviewer Name – Enter the name of the person who Reviewed\QC’ed the
- 1034 data into cell A32
- 1035 ▪ Date – Enter the date on which the data was reviewed into cell D32

1036 11.5.3 Determination of Outliers

1037 The Study Director will use good statistical judgment for determining “unusable” wells that will
 1038 be excluded from the data analysis and will provide an explanation in the study notebook for any
 1039 excluded data. This judgment for data acceptance will include Q-test analysis.

1040 The formula for the Q test is:

1041
$$\frac{\text{Outlier} - \text{Nearest Neighbor}}{\text{Range (Highest} - \text{Lowest)}}$$

1042 where the outlier is the value proposed for exclusion, the nearest neighbor is the value closest to
 1043 the outlier, and the range is the range of the three values (Q values for samples sizes from 3 to 10
 1044 are provided in **Table 11-2**). For example, if the value of this ratio is greater than 0.94 (the Q
 1045 value for the 90% confidence interval for a sample size of three) or 0.76 (the Q value for the 90%
 1046 confidence interval for a sample size of four), the outlier may be excluded from data analysis.

1047 **Table 11-2 Q Test Values**

Number Of Observations	Q Value
2	-
3	0.94
4	0.76
5	0.64
6	0.56
7	0.51
8	0.47
9	0.44
10	0.41

1048

1049 For E2 reference standard replicates (sample size of two), any adjusted RLU value for a replicate
 1050 at a given concentration of E2 is considered and outlier if its value is more than 20% above or
 1051 below the adjusted RLU value for that concentration in the historical database.

1052 11.5.4 Acceptance Criteria1053 11.5.4.1 *Range Finder Testing*

1054 Acceptance or rejection of a test is based on evaluation of reference standard and control results
1055 from each experiment conducted on a 96-well plate. Results are compared to quality controls
1056 (QC) for these parameters derived from the historical database, which are summarized below.

- 1057 • Induction: Plate induction, as measured by dividing the averaged highest E2
1058 reference standard RLU value by the averaged DMSO control RLU value, must
1059 be greater than three-fold.
- 1060 • DMSO control results: Solvent control RLU values must be within 2.5 times the
1061 standard deviation of the historical solvent control mean RLU value.

1062 An experiment that fails either acceptance criterion will be discarded and repeated.

1063 11.5.4.2 *Comprehensive testing*

1064 Acceptance or rejection of a test is based on evaluation of reference standard and control results
1065 from each experiment conducted on a 96-well plate. Results are compared to quality controls
1066 (QC) for these parameters derived from the historical database, which are summarized below.

- 1067 • Induction: Plate induction, as measured by dividing the averaged highest E2
1068 reference standard RLU value by the averaged DMSO control RLU value, must
1069 be greater than three-fold.
- 1070 • Reference standard results: The E2 reference standard concentration-response
1071 curve should be sigmoidal in shape and have at least three values within the linear
1072 portion of the concentration-response curve.
- 1073 • DMSO control results: DMSO control RLU values must be within 2.5 times the
1074 standard deviation of the historical solvent control mean RLU value.
- 1075 • Positive control results: Methoxychlor control RLU values must be above the line
1076 representing the DMSO mean plus three times the standard deviation from the
1077 DMSO mean.


1078 An experiment that fails any single acceptance criterion will be discarded and repeated.

1079

1080 **12.0 RANGE FINDER TESTING**

1081 Agonist range finding for coded substances consists of seven point, 1:10 serial dilutions, with
 1082 each concentration tested in duplicate wells of the 96-well plate. **Figure 12-1** contains a template
 1083 for the plate layout to be used in agonist range finder testing.

1084 **Figure 12-1 Agonist Range Finder Test Plate Layout**



- **Four Point E2 Reference Standard**
- **DMSO (Solvent Control)**
- **Range Finder for Sample #1**
- **Range Finder for Sample #2**
- **Range Finder for Sample #3**
- **Range Finder for Sample #4**
- **Range Finder for Sample #5**
- **Range Finder for Sample #6**

1085


1086 Evaluate whether range finder experiments have met the acceptance criteria
 1087 (see **Section 11.5.4.1**).

1088 To determine starting concentrations for comprehensive testing use the following criteria:

- 1089
- 1090
- 1091
- 1092
- 1093
- 1094
- 1095
- 1096
- 1097
- 1098
- 1099
- 1100
- 1101
- 1102
- 1103
- 1104
- 1105
- 1106
- 1107
- 1108
- 1109
- 1110
- 1111
- 1112
- 1113
- 1114
- 1115
- 1116
- 1117
- If results in the range finder test suggest that the test substance is negative for agonist activity (i.e., if there are no points on the test substance concentration curve that are greater than the mean plus three times the standard deviation of the DMSO control, see **Figure 12-2**), comprehensive testing will be conducted using an 11 point 1:2 serial dilution starting at the maximum soluble concentration.
 - If results in the range finder test suggest that the test substance is negative for agonist activity (i.e., if there are no points on the test substance concentration curve that are greater than the mean plus three times the standard deviation of the DMSO control), and the higher concentrations in the range finder are cytotoxic, comprehensive testing will be conducted using an 11 point 1:2 serial dilution with the lowest cytotoxic concentration as the starting concentration (see **Figure 12-3**).
 - If results in the range finder test suggest that the test substance is positive for agonist activity (i.e., if there are points on the test substance concentration curve that are greater than the mean plus three times the standard deviation of the DMSO control), the starting concentration to be used for the 11-point dilution scheme in comprehensive testing should be one log higher than the concentration giving the highest adjusted RLU value in the range finder. The 11-point dilution scheme will be based on either 1:2 or 1:5 dilutions according to the following criteria:
 - An 11-point 1:2 serial dilution should be used if the resulting concentration range (note: an 11-point 1:2 serial dilution will cover a range of concentrations over approximately three orders of magnitude [three logs]) will encompass the full range of responses based on the concentration response curve generated in the range finder test (see **Figure 12-4**).
 - If the concentration range that would be generated with the 1:2 serial dilution will not encompass the full range of responses based on the concentration response curve in the range finder test (see **Figure 12-5** and **12-6**), an 11-point 1:5 serial dilution should be used instead.

- 1118 • If a substance exhibits a biphasic concentration response curve in the range finder
- 1119 test, both phases should also be resolved in comprehensive testing. In order to
- 1120 resolve both curves, the starting concentration should be based on the peak
- 1121 associated with the higher concentration and should be one log higher than the
- 1122 concentration giving the highest adjusted RLU value in the range finder. As an
- 1123 example, an 11-point 1:5 serial dilution should be used based on the range finder
- 1124 results presented in **Figure12-7**.

1125 **Figure 12-2 Agonist Range Finder (example 1)**




1126 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.

1128

1129


1129 **Figure 12-3 Agonist Range Finder (example 2)**



1130

1131 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.


1132 **Figure 12-4 Agonist Range Finder (example 3)**



1133

1134 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.


1135 **Figure 12-5 Agonist Range Finder (example 4)**



1136

1137 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.

1138 **Figure 12-6 Agonist Range Finder (example 5)**




1139

1140 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.

1141

1141 **Figure 12-7 Agonist Range Finder (example 6)**



1142
 1143 The solid horizontal line represents the mean plus three times the standard deviation of the DMSO control.
 1144

1145 **13.0 COMPREHENSIVE TESTING**

1146 Agonist comprehensive testing for coded substances consists of 11 point serial dilutions (either
 1147 1:2 or 1:5 serial dilutions based on the starting concentration for comprehensive testing criteria in
 1148 **Section 12.0**) with each concentration tested in triplicate wells of the 96-well plate. **Figure 13-1**
 1149 contains a template for the plate layout to be used in agonist comprehensive testing.

1150

1150 **Figure 13-1 Agonist Comprehensive Test Plate Layout**

1151

1152

1153

1154

1155

1156

1157


1158

1159

1160

1161

1162



1163

 **11 Point Duplicate E2 Reference Standard**

1164

 **DMSO (Solvent Control)**

1165

 **Test Substance #1**

1166

 **Test Substance #2**

1167

 **Methoxychlor Control**

1168

1169

1170 Evaluate whether comprehensive experiments have met acceptance criteria (see **Section 11.6.4**)

1171 and graph the data as described in the NICEATM Prism® users guide.

1172 • If the substance has been tested up to the limit dose or the maximum soluble dose,
1173 without causing a significant decrease in cell viability, and there are no points on
1174 the concentration curve that are greater than the mean plus three times the
1175 standard deviation of the DMSO control, the substance is considered negative for
1176 agonism

1177 • If the substance has a positive response (See **Section 6.0**) at any concentration,
1178 the substance is considered positive for agonism.

1179 14.0 COMPILATION OF THE HISTORICAL QUALITY CONTROL DATABASE

1180 Historical databases are maintained in order to ensure that the assay is functioning properly.
1181 Historical databases are compiled using Excel® spreadsheets and are separate from the
1182 spreadsheets used to collect the data for individual test plates. Reference standard and control
1183 data are used to develop and maintain the historical database and are used as quality controls to
1184 determine acceptance of individual test plates.

1185 The sources of the data needed to compile the historical database for the DMSO control are the
1186 experiment specific Excel® data collection and analysis spreadsheets used for LUMI-CELL® ER
1187 agonist and antagonist testing (see **Section 11.5.2** of the LUMI-CELL® ER agonist protocol and
1188 **Section 13.5.2** in the LUMI-CELL® ER antagonist protocol).

1189 14.1 DMSO Control

1190 Open the combined agonist and antagonist LUMI-CELL® ER historical database Excel®
1191 spreadsheet (LUMI_AgandAntQC.xls) and save under a new name using the Excel® “Save As”
1192 function, adding the laboratory designator to the file name (e.g., for Laboratory H, the new name
1193 would be HLUMI_AgandAntQC.xls). Enter the date and experiment name into worksheet
1194 columns A and B respectively. Enter the experimental mean DMSO control value (from cell H37
1195 in the RAW DATA worksheet of the agonist and antagonist Excel® data collection and analysis
1196 spreadsheet) into worksheet column C. Acceptance or rejection of the plate DMSO control data
1197 for range finding and comprehensive testing is based on whether the mean plate DMSO RLU
1198 value falls within 2.5 times the standard deviation of the DMSO value in the historical database
1199 (columns G and H in the DMSO worksheet).

1200

1201 15.0 QUALITY TESTING OF MATERIALS

1202 All information pertaining to the preparation and testing of media, media supplements, and other
1203 materials should be recorded in the Study Notebook.

1204

1204 **15.1 Tissue Culture Media**

1205 Each lot of tissue culture medium must be tested in a single growth flask of cells before use in
1206 ongoing tissue culture or experimentation (**note:** each bottle within a given lot of
1207 Charcoal/Dextran treated FBS must be tested separately).

- 1208 1. Every new lot of media (RPMI and DMEM) and media components (FBS,
1209 Charcoal/Dextran treated FBS, and L-glutamine) must first be tested on the
1210 LUMI-CELL® ER assay prior to being used in any GLP acceptable assays.
- 1211 2. Add 4 µL of DMSO (previously tested) into four separate 13 mm tubes.
- 1212 3. Add 400 µL media (to be tested) to the same tubes.
- 1213 4. Dose an experimental plate as in **Section 12.0**, treating the media being tested as a
1214 test substance.
- 1215 5. Analyze 96-well plate as described in **Section 12.0**, comparing the data from the
1216 DMSO controls made using previously tested tissue culture media to the new
1217 media being tested.
- 1218 6. Use the agonist historical database to determine if the new media with DMSO lies
1219 within 2.5 standard deviations of the mean for the media. If the RLU values for
1220 the new media with DMSO lie within 2.5 standard deviation of the mean for the
1221 historical data on DMSO, the new lot of media is acceptable. If the RLU values
1222 for the new media with DMSO do not lie within 2.5 standard deviations of the
1223 DMSO mean from historical database, the new lot may not be used in the assay.
- 1224 7. Note date and lot number in study notebook.
- 1225 8. If the new bottle passes quality testing as described in **Section 15.1 step 6**, apply
1226 the media to a single flask of cells and observe cell growth and morphology over
1227 the following 2 – 3 days. If there is no change in growth or morphology, the new
1228 media is acceptable for use.

1229

1229 **15.2 G418:**

- 1230 1. New lots of G418 must first be tested on the LUMI-CELL® ER assay prior to
1231 being used in any GLP acceptable assays.
- 1232 2. Add 220 µL of G418 (previously tested) to a single flask containing cells growing
1233 in RPMI.
- 1234 3. Add 220 µL of G418 (to be tested) to a different flask containing cells growing in
1235 RPMI.
- 1236 4. Observe cellular growth and morphology in both tissue culture flasks over a 48 to
1237 72 hour period. If there are no differences in observed growth rate and
1238 morphology between the two flasks, the new G418 lot is acceptable.
- 1239 5. If cellular growth is decreased, or the cells exhibit abnormal morphology, the new
1240 lot of G418 is not acceptable.
- 1241 6. Note date and lot number in study book.

1242 **15.3 DMSO**

- 1243 1. Every new bottle of DMSO must be tested on the LUMI-CELL® ER assay prior
1244 to use in any GLP acceptable assays.
- 1245 2. Add 4 µL of DMSO (to be tested) into four separate 13 mm tubes.
- 1246 3. Add 400 µL media (previously tested) to the same tubes.
- 1247 4. Dose an experimental plate as in **Section 12.0**, treating the DMSO containing
1248 media being tested as a test substance.
- 1249 5. Analyze 96-well plate as described in **Section 12.0**, comparing the data from the
1250 DMSO controls made using previously tested tissue culture media.
- 1251 6. Use the agonist historical database to determine if media with new DMSO lies
1252 within 2.5 standard deviations of the DMSO mean from historical database. If the
1253 RLU values for the media with new DMSO lie within 2.5 standard deviations of
1254 the DMSO mean from the historical database, the new lot of DMSO is acceptable.
1255 If the RLU values for media with new DMSO do not lie within 2.5 standard

1256 deviations of the DMSO mean from historical database, the new lot may not be
1257 used in the assay.

1258 7. Note the date, lot number, and bottle number in study book.

1259 8. If no DMSO has been previously tested, test several bottles as described in
1260 **Section 15.3**, and determine whether any of the bottles of DMSO have a lower
1261 average RLU than the other bottle(s) tested. Use the DMSO with the lowest
1262 average RLU for official experiments.

1263 **15.4 Plastic Tissue Culture Materials**

1264 1. Grow one set of cells, plate them for experiments on plastic ware from the new lot
1265 and one set of cells in the plastic ware from a previous lot, and dose them with E2
1266 reference standard and controls.

1267 2. Perform the LUMI-CELL® ER experiment with both sets of cells.

1268 3. If all of the analysis falls within acceptable QC criteria, then the new
1269 manufacturer's products may be used.

1270

1270 **16.0 REFERENCES**

- 1271 Eli Lilly and Company and National Institutes of Health Chemical Genomics Center. 2005.
1272 Assay Guidance Manual Version 4.1. Bethesda, MD: National Institutes of Health. Available:
1273 http://www.ncgc.nih.gov/guidance/manual_toc.html [accessed 05 September 2006]
- 1274 ICCVAM. 2001. Guidance Document on Using *In Vitro* Data to Estimate *In Vivo* Starting Doses
1275 for Acute Toxicity. NIH Pub. No. 01-4500. Research Triangle Park, NC: National Institute of
1276 Environmental Health Sciences. Available: [http://iccvam.niehs.nih.gov/methods/invidocs/
1277 guidance/iv_guide.pdf](http://iccvam.niehs.nih.gov/methods/invidocs/guidance/iv_guide.pdf) [accessed 31 August 2006]