

Appendix B1

Protocol for CV1 + hAR + Luciferase Assay

**(Provided by Dr. Elizabeth M. Wilson, Departments of
Pediatrics and of Biochemistry and Biophysics,
University of North Carolina, Chapel Hill, NC, USA)**

[This page intentionally left blank]

6 cm dish CV1 LUCIFERASE ASSAY (hAR)

Revised 1-25-02

Monday

1. recount cells: best to count 50-100 cells per 5x5 in hymocytometer, count 2 5x5 grids and average, count $\times 10^4$ = cells/ml
2. plate 0.42×10^6 CV1 cells/6 cm dish by preparing large mix of cells and media so 4 ml media/plate containing 5% bovine calf serum, DMEM-H/20 mM Hepes (2 M Hepes stock, pH 7.2, filter), penicillin and streptomycin, 2 mM L-glutamine, spread cells evenly. Cells usually added from master mix; swirl often while adding cells to the plates.

Tuesday

prepare CaPO_4 precipitates for groups of up to 6 plates using freshly made solutions: for large assays of same DNA, pool the precipitates before adding to the plates.

(a) make 2 M CaCl_2 : 2.94 g $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ bring to 10 ml with sterile autoclaved dH_2O , filter sterilize

(b) make 2X HBS: 500 ml 8.2 g NaCl
 12.5 ml 2 M Hepes
 0.2 g $\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$

bring to 500 ml with sterile autoclaved dH_2O from TC room, pH with 5 N NaOH (takes 150-200 μl), pH to 7.11-7.14, sterile filter, make 27.5 ml aliquots (25 ml needed for 100 6 cm plates), store frozen at -20°C

[for 50 ml 2XHBS: 14 ml 1 M NaCl, 0.25 ml 2 M Hepes Na salt, 750 μl 0.1 M Na_2HPO_4 , bring to 45 ml with sterile ddH_2O (use sterile autoclaved water), add about 45 μl 5 M NaOH, pH to 7.11-7.14, filter sterilize, store pH electrode in pH 7 buffer, NOT H_2O]

(c) prepare DNA one or more days before assay

make dilutions of DNA stocks so additions are ~ 3 -10 μl

add expression and reporter DNA to bottom of 14 ml polystyrene round-bottom (17x100 mm) Falcon tubes, store frozen -20°C

50 ng pCMVhAR (or 10 ng pCMVhAR1-660 (ABC))

5 μg MMTV-luciferase

[For PSA-luciferase use 5 μg reporter/dish, 100 ng pCMVhAR/dish]

(d) per 6 cm dish, add to tubes containing DNA:

210 μl sterile H_2O

30 μl 2 M CaCl_2 (final 0.125 M CaCl_2)

240 μl 2X HBS, vortex briefly, let sit 30 min at RT

vortex briefly, add 475 μl of mix per well, return plates to incubator, incubate 4 h

example for 6 dishes:

<u>DNA</u>	<u>H_2O</u>	<u>2 M CaCl_2</u>	<u>2XHBS</u>
0.3 μg pCMVhAR	1.26 ml (2 x 630 μl)	180 μl	1.44 ml (2 x 720 μl)
30 μg MMTV-Luc			

(e) aspirate plates, add 1.5 ml glycerol shock medium, incubate 3 min RT, aspirate, wash 4 ml PBS, aspirate, add 4 ml serum free, phenol red free DMEM-H, Hepes, P/S, glutamine \pm hormone, return to incubator for overnight

Glycerol Shock Medium: use 5% DMEM-H red

dishes \times 1.5 ml/dish = total volume (make extra)

total volume \times 15% = amount of glycerol

total volume – amount of glycerol = amount of media to add with glycerol

Wednesday aspirate media, add 4 ml fresh phenol red-free, serum-free DMEM-H, P/S, Hepes, glutamine \pm hormone, add DHT to stock media as needed and add to plates, incubate overnight 37°C

Thursday remove media, wash with 4 ml PBS, aspirate to dry; add 0.5 ml/plate lysis buffer; rock plates 20-30 min, not much longer at RT. Transfer 100 μ l from each well to 96 well Nunc flat bottom standard nontreated white plate. Read on LumiStar automated luminometer that injects 100 μ l luciferin stock and 100 μ l reading buffer

Lysis buffer: 2 mM EDTA, 1% Triton X-100, 25 mM Trizma (Tris base) phosphate, pH 7.8

D-Luciferin: prepare 1 mM D-luciferin (K⁺ salt, MW 318.41) in dH₂O store in 10 ml aliquots at -20°C covered with foil (D-luciferin is light sensitive) use 100 μ l/sample, save extra at -20°C, Na⁺ salt sometimes turns yellow but is probably still good, pH of H₂O might be off, better to use K⁺ salt) (from Analytical Luminescence)

Reading buffer: The optimal pH for the reaction is pH 7.8 cold; if glycylglycine and ATP are carefully pH cold, then the final will be pH 7.8

Stock	Amount to <u>20 ml final</u>	Amount to <u>100 ml final</u>	<u>Final conc</u>
0.5 M glycylglycine, pH 7.8 cold	1 ml	5 ml	25 mM
1 M MgCl ₂	300 μ l	1.5 ml	15 mM
100 mM ATP in dH ₂ O	1 ml	5 ml	5 mM
(bring to pH 7.8 cold with 1 M NaOH, CRITICAL, store -80°C 1 ml aliquots)			
Sigma - tissue culture grade			
50 mg/ml BSA dH ₂ O	200 μ l	1 ml	0.5 mg/ml
dH ₂ O	17.5 ml	87.5 ml	
need 40 ml for 100 plates			