

# **Blood levels of folate over time, current US levels, and differences between assessment methods**

**Christine M Pfeiffer, PhD**

Chief, Nutritional Biomarkers Branch

NTP/NIEHS Expert Panel Meeting “Identifying Research Needs for  
Assessing Safe Use of High Intakes of Folic Acid”

May 11-12, 2015 in Bethesda, MD

National Center for Environmental Health  
Division of Laboratory Sciences



## Disclosure

Nothing to disclose

National Center for Environmental Health  
Division of Laboratory Sciences




## **Presentation outline**

- **Assays for serum and red blood cell (RBC) folate**
  - Brief overview of available methods
  - Issues with (lack of) comparability of data
  - Issues with cutoff values
- **Serum and RBC folate status pre- vs. post-fortification**
- **Current US blood folate levels**
  - Post-fortification concentrations of serum and RBC folate
  - Post-fortification prevalence estimates of low blood concentrations
  - Post-fortification concentrations of serum folate forms, including unmetabolized folic acid (UMFA)
  - Factors associated with higher UMFA concentrations



## ASSAYS FOR SERUM AND RBC FOLATE



## Folate structure allows variations in three areas



## Serum and RBC folate – 2012 UK NEQAS PT program



## Main laboratory methods for serum and RBC folate

Method type	Advantages	Disadvantages
<u>Microbiologic assay (MBA)</u> for total folate	<ul style="list-style-type: none"> <li>• Small sample volume</li> <li>• Inexpensive</li> <li>• Suited for low resource setting</li> <li>• Measures all biologically active forms approximately equally (however, calibration with 5-methylTHF generates lower results than calibration with folic acid)</li> </ul>	<ul style="list-style-type: none"> <li>• Relatively laborious and takes 2 days to result</li> <li>• Replicates needed due to higher imprecision</li> <li>• Multiple dilutions needed due to limited linear range</li> <li>• Inhibited by presence of antibiotics or antifolates</li> </ul>
<u>Competitive protein binding assay (CPBA)</u> for total folate	<ul style="list-style-type: none"> <li>• User friendly and minimum operator involvement</li> <li>• High sample throughput</li> <li>• Suited for clinical setting</li> <li>• Generally good precision (~5%)</li> </ul>	<ul style="list-style-type: none"> <li>• Questionable accuracy due to different affinities of folate forms to FBP</li> <li>• Less suited for long-term studies due to potential lot-to-lot variability</li> <li>• Matrix effects likely with sample dilution</li> </ul>
<u>Chromatography-based assay</u> for folate (various detectors; recently MS/MS)	<ul style="list-style-type: none"> <li>• Information on folate vitamers</li> <li>• Highly selective and specific</li> <li>• Good analytical sensitivity and precision</li> <li>• Suited for research setting</li> </ul>	<ul style="list-style-type: none"> <li>• Expensive instrumentation and technical service, experienced operator</li> <li>• Complex sample extraction/clean-up</li> <li>• Conversion of polyglutamates to monoglutamates needed for whole blood</li> <li>• Summation of folate forms to total folate</li> </ul>



## Folate methods and data in NHANES

NHANES	Serum folate	RBC folate
<b>Pre-fortification</b>		
1988-1994	BioRad RIA	BioRad RIA
<b>Post-fortification</b>		
1999-2000	BioRad RIA	BioRad RIA
2001-2002	BioRad RIA	BioRad RIA
2003-2004	BioRad RIA	BioRad RIA
2005-2006	BioRad RIA	BioRad RIA
2007-2008	MBA (LC-MS/MS)	MBA
2009-2010	MBA	MBA
2011-2012	LC-MS/MS	MBA
2013-2014	LC-MS/MS	MBA
2015-2016	LC-MS/MS	MBA

Converted to MBA-equivalent data  
*Pfeiffer et al. J Nutr 2012*


MBA approx equivalent to LC-MS/MS  
for serum folate  
*Fazili et al. Clin Chem 2007*  
*Yetley et al. AJCN 2011*




## Folate methods and data in the literature

- MBA, CPBA, chromatography-based – little information on how the various methods compare at any given point in time
- MBA:
  - Limited number of labs using this assay
  - Assay results may vary depending on microorganism (antibiotic-resistant vs. wild type) and calibrator (folic acid vs. 5-methylTHF vs. 5-formylTHF)
  - Most comprehensive comparison data: *Pfeiffer et al. J Nutr 2011*
  - 2015 CDC Folate Round Robin for microbiologic assay labs
- CPBA:
  - Comparison data from a few studies, but question whether assay changed over time (*Gunter et al. Clin Chem 1996; Pfeiffer et al. Clin Chem 2001; Owen et al. Am J Clin Pathol 2003; Clifford et al. J Nutr 2005*)
- LC-MS/MS:
  - No comparison data available
  - 2015 CDC Serum Folate Round Robin for LC-MS/MS labs

## Cutoff values for “abnormal” folate levels





## SERUM AND RBC FOLATE STATUS PRE- VS. POST-FORTIFICATION



\*BioRad data


CDC's 2012 Second National Nutrition Report at [www.cdc.gov/nutritionreport](http://www.cdc.gov/nutritionreport)

**Post-fortification serum and RBC folate concentrations were  
~2.5x and 1.5x pre-fortification concentrations, respectively**





\* MBA-equivalent data

### Higher serum and RBC folate 95<sup>th</sup> percentile concentrations for supplement users vs. non-users



\* MBA-equivalent data

## Relationship between blood folate concentrations and total vitamin intake quartiles in adult supplement users and non-users NHANES 2003-2006



\*BioRad data



# CURRENT BLOOD FOLATE LEVELS



*British Journal of Nutrition*, page 1 of 13  
doi:10.1017/S0007114515001142  
© The Authors 2015. This is a work of the U.S. Government and is not subject to copyright protection in the United States

## Folate status and concentrations of serum folate forms in the US population: National Health and Nutrition Examination Survey 2011–2

Christine M. Pfeiffer<sup>1\*</sup>, Maya R. Sternberg<sup>1</sup>, Zia Fazili<sup>1</sup>, David A. Lacher<sup>2</sup>, Mindy Zhang<sup>1</sup>, Clifford L. Johnson<sup>2</sup>, Heather C. Hamner<sup>3</sup>, Regan L. Bailey<sup>4</sup>, Jeanne I. Rader<sup>5</sup>, Sedigheh Yamini<sup>5</sup>, R. J. Berry<sup>3</sup> and Elizabeth A. Yetley<sup>1</sup>

<sup>1</sup>National Center for Environmental Health, Centers for Disease Control and Prevention, Atlanta, GA, USA

<sup>2</sup>National Center for Health Statistics, Centers for Disease Control and Prevention, Hyattsville, MD, USA

<sup>3</sup>National Center for Birth Defects and Developmental Disabilities, Centers for Disease Control and Prevention, Atlanta, GA, USA

<sup>4</sup>Office of Dietary Supplements, National Institutes of Health, Bethesda, MD, USA

<sup>5</sup>Center for Food Safety and Applied Nutrition, Food and Drug Administration, College Park, MD, USA

## U.S. Women of Childbearing Age Who Are at Possible Increased Risk of a Neural Tube Defect-Affected Pregnancy Due to Suboptimal Red Blood Cell Folate Concentrations, National Health and Nutrition Examination Survey 2007 to 2012

Sarah C. Tinker<sup>\*1</sup>, Heather C. Hamner<sup>2</sup>, Yan Ping Qi<sup>1</sup>, and Krista S. Crider<sup>1</sup>



## Current folate status in the US population NHANES 2007 – 2012

Survey period	Serum folate (nmol/L)	RBC folate (nmol/L)
2007 – 2008	39.5 (37.7 – 41.3)	1120 (1070 – 1160)
2009 – 2010	38.2 (37.2 – 39.3)	1040 (1010 – 1070)
2011 – 2012	41.4 (40.1 – 42.9)	1050 (1010 – 1090)

WHO cutoff:			
	10 nmol/L	340 nmol/L	906 nmol/L
Survey period	Serum folate <13.7 nmol/L	RBC folate <624 nmol/L	RBC folate <748 nmol/L
2007 – 2008	3.5%	7.6%	22.8% (women 12-49 y)
2009 – 2010	3.9%	9.4%	
2011 – 2012	1.1%	9.0%	


*Pfeiffer et al. Br J Nutr 2015*

*Tinker et al. 2015*


## **Current blood folate concentrations NHANES 2011-2012**

- NHANES 2011-2012 provides the first data on serum folate forms for persons 1 y and older by demographic and selected physiologic and lifestyle variables
- Concentrations of 5-methylTHF (100%), UMFA (99.9%), MeFox (98.8%), and THF (85.2%) mostly detectable
- 5-FormylTHF (3.6%) and 5,10-methenylTHF (4.4%) rarely detected
- Contribution to total folate: 5-methylTHF (86.7%), UMFA (4%), non-methylfolate (4.7%), MeFox (4.5%)
- All biomarkers showed higher concentrations with recent folic acid-containing supplement use

## SERUM UNMETABOLIZED FOLIC ACID CONCENTRATIONS



### Serum UMFA contributed 4% to total folate, NHANES 2011-2012



Concentrations varied with age, sex, race-ethnicity, fasting status, eGFR, BMI, BSA, serum cotinine, alcohol intake, and folic acid supplement use

The Journal of Nutrition. First published ahead of print December 10, 2014 as doi: 10.3945/jn.114.201210.

The Journal of Nutrition  
Nutritional Epidemiology



## **Unmetabolized Folic Acid Is Detected in Nearly All Serum Samples from US Children, Adolescents, and Adults<sup>1-4</sup>**

Christine M Pfeiffer,<sup>5\*</sup> Maya R Sternberg,<sup>5</sup> Zia Fazili,<sup>5</sup> Elizabeth A Yetley,<sup>6</sup> David A Lacher,<sup>7</sup> Regan L Bailey,<sup>6</sup> and Clifford L Johnson<sup>7</sup>

<sup>5</sup>National Center for Environmental Health, CDC, Atlanta, GA; <sup>6</sup>Office of Dietary Supplements, NIH, Bethesda, MD; and <sup>7</sup>National Center for Health Statistics, CDC, Hyattsville, MD

- Prevalence of UMFA >1 nmol/L was 33% overall and 21% in fasting (≥8 h) adults in NHANES 2007-2008
- UMFA >1 nmol/L was largely but not entirely explained by fasting status and by total folic acid intake from diet and supplements

## Summary

- Folate assays have not yet been standardized and results are not comparable across assays or laboratories
- Fortification has significantly increased blood folate levels in the US population
- Post-fortification prevalence of deficient blood folate levels is <10%
- Post-fortification blood folate levels have been fairly constant over a period of ~15 years
- Post-fortification detection of serum UMFA is nearly ubiquitous and concentrations >1 nmol/L are largely but not entirely explained by fasting status and by total folic acid intake from diet and supplements





## **Acknowledgments**

Staff of the Nutritional Biomarkers Branch  
Colleagues at NCHS/NHANES, NIH/ODS, and FDA



## Questions?

**For more information please contact Centers for Disease Control and Prevention**

1600 Clifton Road NE, Atlanta, GA 30333

Telephone: 1-800-CDC-INFO (232-4636)/TTY: 1-888-232-6348

Visit: [www.cdc.gov](http://www.cdc.gov) | Contact CDC at: 1-800-CDC-INFO or [www.cdc.gov/info](http://www.cdc.gov/info)

The findings and conclusions in this report are those of the authors and do not necessarily represent the official position of the Centers for Disease Control and Prevention.

National Center for HIV/AIDS, Viral Hepatitis, STD, and TB Prevention  
Division of Laboratory Sciences

