

ToxCast / Tox21 Data Overview

Richard Judson EPA National Center for Computational Toxicology

COMPUTATIONAL TOXICOLOGY

The views expressed in this presentation are those of the author and do not necessarily reflect the views or policies of the U.S. Environmental Protection Agency.



A Problem Driving CompTox Development

- ~1,000 chemicals most people are often exposed to
- ~10,000 chemicals many people are exposed to at least occasionally
- ~100,000 chemicals some people are exposed to at least occasionally
- Some of these chemicals / exposures are responsible for
 - -Cancers
 - -Birth defects
 - Reproductive disorders
 - Both human and ecological species
- Can we prioritize chemicals for further testing without knowing everything?




ToxCast Projects

- ToxCast / Tox21 is a large-scale in vitro screening program
 - -ToxCast: 1,000 chemicals in ~500 assays
 - -Tox21: 10,000 chemicals in ~50 assays (drugs+EPA+NTP chemicals)
- Reverse toxicokinetics (RTK) approach for estimating maximum "no bioactivity" dose estimate
- Input for Endocrine Screening Prioritization and other programs
- Models for Prioritization of Targeted Testing
 - -Test examples have known MOA, activity
 - Goal is to qualify in vitro assays for use in prioritizing testing of untested compounds



United States Environmental Protection The ToxCast Phase I Chemicals (320)

309 Unique Structures Replicates for QC





ToxCast Assays

~500 Total Endpoints

Cellular Assays

Biochemical Assays

- Protein families
 - GPCR
 - NR
 - Kinase
 - Phosphatase
 - Protease
 - Other enzyme
 - Ion channel
 - Transporter
- Assay formats
 - Radioligand binding
 - Enzyme activity
 - Co-activator recruitment

Primarily Human / Rat

Cell lines

- HepG2 human hepatoblastoma
- A549 human lung carcinoma
- HEK 293 human embryonic kidney

Primary cells

- Human endothelial cells
- Human monocytes
- Human keratinocytes
- Human fibroblasts
- Human proximal tubule kidney cells
- Human small airway epithelial cells
- Rat hepatocytes
- Mouse embryonic stem cells (Sid Hunter)

Biotransformation competent cells

- Primary rat hepatocytes
- Primary human hepatocytes

Assay formats


- Cytotoxicity
- Reporter gene
- Gene expression
- Biomarker production
- High-content imaging for cellular phenotype




Data Analysis: What is a hit?




Attagene ERa_TRANS







The ToxCast In Vitro Data Set



Many hits – median value=50


Fewer in cell-free HTS (Novascreen, red)

Many hits are at or near top of tested concentration range

A few are in nanomolar range




Reference Toxicity Database (ToxRefDB)






Chemicals show a wide range of activity

Minimal Pathways, 30 uM







Where is the data?

- ToxCastDB contains the ToxCast and ToxRefDB data for browsing
 - -http://actor.epa.gov/actor/faces/ToxMiner/Home.jsp
- ACToR contains a much large amount of data on other types of toxicity and exposure
 - http://actor.epa.gov




ACTOR

Aggregated Computational Toxicology Resource

http://actor.epa.gov/







Data Definitions

Substance

- A chemical from one source
- Name(s), CASRN
- Source-specific unique ID
- -Assay Data

Compound

- -Chemical structure from one source
- -Source-specific unique ID

Generic Chemical

- -CASRN
- Link to many substance (each with same CASRN)
- Link to at most one compound
- -Links to all assay data from substances with same CASRN



Data Definitions

Assay

- A collection of data on one or more substances
- -Comes from one data source
- Can have several types of data included
- Looks like an Excel spreadsheet

Assay Component

One column of an assay table

Assay Result

A data value for one substance and one assay component



Data Definitions

- Assay Phenotype
 - Type of disease associated with the assay
 - Carcinogenicity, GeneTox, ...
- Assay Category
 - -Type of data: tabular, links to the web, human exposure
 - Allows assays to be grouped together
- Data Collection
 - A source of data
 - -Substances
 - Compounds
 - –Assays




Data Sets

- 1. Attagene
- 2. ACEA
- 3. BioSeek
- 4. Cellumen
- 5. Gentronix
- 6. NCGC
- 7. Novascreen
- 8. CellzDirect
- 9. Solidus
- 10.ToxRefDB

- Perturbation Scores
 - -Genes
 - -Pathways
 - Ingenuity
 - KEGG
 - PathwayCommons
- PhysChem Properties
 - -EPISuite
 - -LeadScope
 - -QikProp
- Structure Classifiers




Multiple Assays per Endpoint





Browsing ToxCast / Tox21 Data


http://actor.epa.gov/actor/faces/ToxMiner/Home.jsp



Home page – search by name or CASRN



ToxCastDB Chemical Page






ToxCastDB Assay Page

82657-04-3 37.0


80-05-7



Bifenthrin Risphenol A




ToxCastDB Data Collection Page






Genes associated with assays

Lists all assays associated with each gene NOTE: this page is slow to load





ACToR Chemical Page – many sources of data





Coming Soon

- Mouse Embryonic Stem Cell Data (With Sid Hunter, EPA NHEERL)
- RTK Data Reverse Toxicokinetics (with R. Thomas, Hamner)
- ExpoCastDB (With collaborators in EPA NERL)



Outline: ToxCast Data Sets Overview

- Data Sets: Bio-centric View
 - Assay overview
 - -How were the data calculated?
 - –How did the replicates perform?
 - Example results
 - –What did the assay measure?
 - Potential assay artifacts



ToxCast Phase I Data Sets

- Biochemical
 - -Novascreen

- Cell-Based
 - -ACEA
 - Attagene
 - -BioSeek
 - -Cellumen
 - -CellzDirect
 - -Gentronix
 - -NCGC
 - -Solidus



Novascreen: 239 Biochemical Assays (Abstract 59)

- Protein super-families
 - -GPCR
 - Kinase
 - Phosphatase
 - Protease
 - Ion channel
 - Nuclear receptor
 - Other enzyme
 - CYP P450 inhibition

- Various formats:
 - Radioligand receptor binding
 - Fluorescent receptor binding
 - Fluorescent enzyme substrateintensity quench
 - Fluorescent enzyme substratemobility shift
- Initial screening:
 - $-25 \mu M$ in duplicate
 - 10 μM in duplicate (CYPs)
- Normalize data to assay window
 - % of control activity (central reference – scalar reference)



Novascreen: What do the assays measure?

- Mainly direct effects of chemical on target protein
 - Enzyme activity
 - Ligand binding
- False positives:
 - Fluorescent compounds—fluorescing and quenching
 - Reactive compounds/covalent modification of target
 - Physical effects—colloid aggregation of target
 - Operational
- False negatives:
 - Solubility
 - Inappropriate assay conditions
 - Operational
 - Target protein not physiological
 - Lack of biotransformation

Data Correction Algorithm Examples (Additive)

		AChE											Caspase 10													
		1	2	3	4	5	6	7	8	9	10	11	12		1	2	3	4	5	6	7	8	9	10	11	12
normalized	А	-3	3	-29	-100	-98	-100	-68	-68	-89	-89	-98	-98	Α	25	0	-3	-97	-100	-101	-35	-33	-76	-80	-100	-99
	В	-4	7	-29	-12	-10	-12	-27	-22	-30	-30	-25	-16	В	9	-17	-16	-23	-28	-17	-13	-16	-22	-15	-22	-3
	С	8	2	-23	-23	-19	-26	-25	-26	-26	-28	-25	-13	С	3	-22	-10	-27	-24	-20	-39	-31	-24	-28	-26	-6
	D	-14	1	-17	-21	-23	-25	-18	-26	-26	-27	-25	-17	D	15	-22	-15	-27	-23	-27	-13	-15	-17	-18	-20	-3
	Ε	-11	-5	-23	-25	-16	-20	-26	-22	-27	-30	-26	-13	Ε	3	-16	-16	-17	-18	-7	-6	-13	2	5	-17	12
	F	-16	0	-7	-7	-17	-21	-22	-24	-29	-33	-19	-17	F	-11	-34	-18	-26	-29	-21	-22	-19	-18	-19	-25	0
	G	-5	4	-32	-18	-25	-17	-25	-28	-25	-32	-20	-13	G	-23	-30	-27	-24	-17	-6	-8	-11	-22	-12	-9	10
	Н	-85	-78	-22	-22	2	-14	-17	-17	NaN	NaN	NaN	NaN	Н	-15	-28	-17	-20	-14	-28	-5	4	-19	-4	-18	5
		1	2	3	4	5	6	7	8	9	10	11	12		1	2	3	4	5	6	7	8	9	10	11	12
Normalized & corrected	А	-1	1	-19	-100	-98	-100	-67	-67	-89	-89	-98	-98	Α	23	0	-3	-97	-100	-101	-36	-34	-76	-81	-100	-99
	В	0	7	-20	-10	-2	0	2	9	-15	-16	-9	-4	В	13	6	7	0	-6	-3	3	2	-6	1	-5	0
	С	10	10	-3	-3	0	-6	-2	-5	-16	-18	-11	-1	С	6	-1	10	-7	-4	-5	-21	-13	-5	-10	-7	-5
	D	4	14	7	2	2	-5	6	0	-13	-11	-9	-3	D	18	0	2	-7	0	-7	6	4	0	-1	-4	1
	Ε	7	13	4	0	7	4	-2	2	-15	-18	-10	-1	Ε	1	1	2	-3	-5	-2	7	0	14	17	-6	9
	F	5	22	16	16	9	6	3	3	-15	-17	-3	-2	F	-9	-12	4	-3	-11	-9	-4	-1	-7	-4	-11	-1
	G	12	23	-4	6	1	5	3	0	-11	-16	-1	4	G	-17	-5	-5	-4	-1	-1	13	6	-2	3	5	8
	Н	-52	-44	7	9	20	8	4	5	NaN	NaN	NaN	NaN	Н	-9	-11	2	0	4	-14	9	15	-4	10	-7	8


United States
Environmental Protection
Agency

Novascreen Concentration-Response


Testing

Retest actives:


- Median absolute deviation (MAD) median lx-xmed
 two MADs or 30% activity
- -8 conc/3-fold serial dilutions
 - 50 μM high conc
 - 25 μM high conc for CYPs
- Normalize to assay window
- Fit % Activity data to 3- or 4parameter Hill function
 - Sometimes had to fix top or bottom of curve
 - Did not extrapolate beyond testing range
 - Manual or automated removal of obvious outliers



Novascreen




Novascreen: Example Curve Fits



ACEA: Real Time Cell Growth Kinetics



- Cytotoxicity with potential mechanistic interpretation
- Human A549 lung carcinoma cell line
 - –ACEA experience with line
 - -Reference compound effects
- Concentration-response testing
 - -8 conc/3-fold serial dilutions
 - Duplicate wells
- Real-time measuremens during exposure (0-72 hr)
- IC50 and LELs calculated






ACEA: What is measured?

- General cytotoxicity in transformed cell line
- False positives
 - -No obvious
 - -Operational
- False negatives
 - Operational error
 - –Solubility
 - Lack of appropriate toxicity targets
 - -Lack of biotransformation




ACEA: Data examples

Replicate Analysis:



Example Plots:







Attagene Multiplexed Transcription Factor Assays


- Modulation of TF activity in human hepatoma HepG2 cells
- Multiplexed reporter gene assay
 - -cis 52 assays (response element driving reporter)
 - trans 29 assays (GAL4-NR_LBD driving reporter) "ligand detection"
- IC50 for cytotoxicity measured first in HepG2
- High concentration either 100 μM or 1/3 calculated IC50 for cytotoxicity
- Seven concentrations, 3-fold serial dilutions, 24 hr exposure
- Cells harvested, RNA isolated, processed for reporter gene quantitation
- LEL provided in data set




Attagene Technology


Cis: AhR

Multiplexed Reporter Gene Assay





■ TV000392





Attagene: What Is Being Measured?

cis Assays


- Up/down regulation of endogenous transcription factor activity in transformed cell line
- False positives
 - General cytotoxic response resulting in non-specific transcriptional activity
 - Promiscuity of response elements
 - Statistical, not biologically, significant response
 - Operational
- False negatives
 - Solubility
 - Cytotoxicity
 - Operational
 - Lack of endogenous machinery
 - Lack of biotransformation


- trans Assays
 - NR agonist activity
 - False positives
 - General cytotoxic response resulting in non-specific transcriptional activity
 - Statistical, not biologically, significant response
 - Operational
 - False negatives
 - Solubility
 - Cytotoxicity
 - Operational
 - · Lack of endogenous machinery
 - Lack of biotransformatioin



Attagene: corresponding *cis* and *trans* assays









cis: ERE

trans: ERa





Office of Research ar National Center for Co


■ TV000003


Attagene: Data calculation challenges


Attagene


United States
Environmental Protection
Agency

- No positive reference compound for each endpoint
- Responses, especially for cis assays, tended to be montonic so potency value difficult to derive
- Biological vs statistical significance for LEL












BioSeek: BioMAP® Technology Platform (Abstract 24)


Assays



Profile Database



Informatics



Human primary cells Disease-like culture conditions Biological responses to drugs and stored in the database

Specialized informatics tools are used to mine and analyze biological data

Primary Human Cell-Based Assay Platform for Human Pharmacology



BioSeek Assays Tested Against ToxCast_320

System		Cell Types	Environment	Readouts
3C	9	Endothelial cells	IL-1β+TNF-α+IFN-γ	MCP-1, VCAM-1, ICAM-1, Thrombomodulin, Tissue Factor, E-selectin, uPAR, IL-8, MIG, HLA-DR, Prolif., Vis., SRB (13)
4H		Endothelial cells	IL-4+histamine	VEGFRII, P-selectin, VCAM-1, uPAR, Eotaxin-3, MCP-1, SRB (7)
LPS		Peripheral Blood Mononuclear Cells + Endothelial cells	TLR4	CD40, VCAM-1,Tissue Factor, MCP-1, E-selectin, IL-1a, IL-8, M-CSF, TNF-a, PGE2, SRB (11)
SAg		Peripheral Blood Mononuclear Cells + Endothelial cells	TCR	MCP-1, CD38, CD40, CD69, E-selectin, IL-8, MIG, PBMC Cytotox., SRB, Proliferation (10)
BE3C		Bronchial epithelial cells	IL-1β+TNF-α+IFN-γ	uPAR, IP-10, MIG, HLA-DR, IL-1a, MMP-1, PAI-1, SRB, TGF-b1, tPA, uPA (11)
HDF3CGF		Fibroblasts	IL-1β+TNF-α+IFN-γ +bFGF+EGF+PDGF-BB	VCAM-1, IP-10, IL-8, MIG, Collagen III, M-CSF, MMP-1, PAI-1, Proliferation, TIMP-1, EGFR, SRB (12)
KF3CT	93	Keratinocytes + Fibroblasts	IL-1β+TNF-α+IFN-γ +TGF-β	MCP-1, ICAM-1, IP-10, IL-1a, MMP-9, TGF-b1, TIMP-2, uPA, SRB (9)
SM3C	3	Vascular smooth muscle cells	IL-1β+TNF-α+IFN-γ	MCP-1, VCAM-1, Thrombomodulin, Tissue Factor, IL-6, LDLR, SAA, uPAR, IL-8, MIG, HLA-DR, M-CSF, Prolif., SRB (14)



BioSeek Assays

- Chemicals tested at 4 concentrations: 40, 13.3, 4.4, 1.5 μM, single well
- Exposure started 1 hr before stimulation of cell signaling pathways
- Following 24 hr exposure, endpoints measured by ELISA (also Alamar blue, SRB staining, and microscopy)

KF3CT

Data calculated vs 6 DMSO controls as fold-change

HDF3CGF

Up and down regulation distinguished


BE3C

Repeats of colchicine positive control

LELs determined

LPS

SAa



Replicate 1

BioSeek



BioSeek Assays: What is Being Measured?


- Effects of chemicals on signaling pathways in primary human cells
- False positives
 - Cytotoxicity (down-regulation in particular)
 - -Statistical vs biological significance
 - -Operational
- False negatives
 - -Solubility
 - Lack of biotransformation
 - -Cytotoxicity
 - Operational

United States Environmental P

Agency

BioSeek MOA Signatures


Defined by reference cmpd behavior in assays used



H₂N O S N S

Myxothiazol (antibiotic) is an inhibitor of coenzyme Q - cytochrome c reductase (aka complex III). It prevents electron transfer from Fe-S center to ubiquinone in the electron transport chain.

Strobilurin fungicides: MOA-blocks electron transfer within respiratory chain




Pyraclostrobin

Trifloxystrobin



Cellumen: High-Content Screening of Cellular Phenotypic Toxicity Parameters (Abstract 38)

- Technology: automated fluorescent microscopy
- Objective: Determine effects of chemicals on toxicity biomarkers in a cell culture of HepG2 and primary rat hepatocytes





Cellumen: CellCiphr™ Cytotoxicity Panel

- 10-point conc-response (200 μM-39 nM) in duplicate
- Three time points (1 hr, 24 hr, 72 hr)
- 11 endpoints per assay

Biomarker	Measurement	Positive Control	Z'
Stress Pathway	Phospho-c-jun	Anisomycin	0.63
Oxidative Stress	Phospho-Histone H2A.X	Camptothecin	0.7
Mitochondrial Function	Mitochondrial membrane potential	CCCP	0.55
Mitochondrial Mass	Mitochondrial mass	CCCP	0.35
Cell Loss	Cell number	Camptothecin	0.56
Cell Cycle	DNA content	Paclitaxel	0.54
DNA Degradation	DNA structure	Paclitaxel	0.6
Nuclear Size	Area of nuclear region	Paclitaxel	0.63
DNA Damage	Detection of p53	Camptothecin	0.43
Mitotic Arrest	Phospho-Histone-H3	Paclitaxel	0.63
Cytoskeletal Integrity	Detection of α -tubulin	Paclitaxel	0.3



Cellumen: What is Being Measured?

- Cellular toxicity phenotypes in a transformed cell line
- False positives
 - Imaging artifacts
 - -Fluorescent compounds
 - -Statistical vs biological significance
 - Operational
- False negatives
 - Solubility
 - -Cytotoxicity
 - Lack of biotransformation
 - -Operational




Cellumen: Data Calculation

- Data not normalized to controls (with a few exceptions)
- Fit to 3- or 4-parameter Hill equation
- AC50 reported with these rules:
 - For Cell Loss, AC50 is reported as calculated
 - For other endpoints, if AC50 for endpoint is > AC50 for Cell Loss at the corresponding exposure time, AC50 for endpoint is set to 100 μ M (to account for imaging artifacts of cytotoxicity

Issues:

- Lack of positive controls for all endpoint/time combinations
- Large differences in maximal response
- noisy curves due in part to effects of cytotoxicity




Cellumen: Data Examples




Cell Loss

Mitochondrial Membrane **Potential**

> **DNA Damage**









CellzDirect: XME Gene Expression in Primary Human Hepatocytes (Abstract 22)

- Primary human hepatocytes from two donors used
- Cells exposed for 6, 24, and 48 hr; medium/chemical refreshed daily
- Concentrations tested: 40, 4, 0.4, 0.04, and 0.004 µM
- 16 Genes measured in multiplexed RNAse protection assay (qNPA)
- Genes targeted XME and transporters



Target Gene Categories
CYP450 (6)
Transporter (4)
Phase II Metabolism (3)
Cholesterol Synthesis (1)
Endogenous control (2)






CellzDirect: What is Being Measured?

- Up/down regulation of mRNA for XME and transporters in primary human hepatocytes
- False positives:
 - -General effect of cytotoxicity on transcriptional activity
 - -Statistical vs biological significance
 - Operational
- False negatives:
 - -Solubility
 - -Cytotoxicity
 - -Operational
 - Lack of biotransformation
 - Inter-individual donor variation




CellzDirect: Data Calculation

- Normalized to solvent control; expressed as fold-change
- Curves fit to Hill equation
 - Upper and lower limits defined by minimum and maximum responses observed over dataset of a particular gene/donor/time
- LELs determined
- Large variations in maximal responses
 - Biological vs statistical significance
- Measuring endogenous promoter activity reflects complex, multifactorial regulation of gene expression
- 6 hr exposure data not provided due to high variability associated with limited time for gene induction




CellzDirect: Data Examples






Gentronix: GADD45a Reporter Gene Assay for DNA Damage (Abstract 41)

- TK6 cell line expressing GFP under control of GADD45a promoter
- Cells exposed at 200, 100, 50 μM for 24 and 48 hr
- Cytotox assay to discount artifacts
- Retested at lower conc if cytotoxic





United States Environmental Protection Agency

Gentronix: Data Calculations


- If induction of GFP fluorescence >50%1 test called genotoxic
- If over 2 or 3 concentrations, strongly genotoxic
- LELs calculated
- Replicate analysis: no actives among replicates

- GADD45a promoter activity in transformed cell line
- False positives:
 - Cytotoxicity resulting in general transcriptional activity
 - Cellular stress other than DNA damage
 - Statistical vs biological validation
 - Operational
- False negatives:
 - Solubility
 - Cytotoxicity
 - Lack of biotransformation
 - Operational




Solidus: Biotransformation Chip and Effect of Cytotoxicity (Abstract 30)


- Alginate-immobilized Phase I and Phase II enzymes
- ToxCast_320 exposed 6 hr to:
 - Control
 - -Ph I
 - -Ph II
 - -Ph I and Ph II
- 960 uM high conc/4-fold serial dilutions/9 concentrations/5 replicates
- Cytotoxicity in Hep3B measured 48 hr later





Solidus

- United States
 Environmental Protection
 Agency
- Data normalized to control values
- Concentration-response data fit to Hill equation
- LC50 determined for each assay condition











Solidus: What is Being Measured?

- Effect of Phase I and Phase II enzymes on cytotoxicity activity of chemicals against a transformed cell line
- False positives:
 - -Possibly not optimized Phase I and/or Phase II mix
 - -Operational
- False negatives
 - -Solubility
 - -Possibly not optimized Phase I and/or Phase II mix
 - Availability of compound from alginate-immobilized enzyme matrix
 - -Operational


United States Environmental Protection Agency

NCGC Reporter Gene Assays

- Nuclear Receptors
 - GAL4 System (ligand detection assay)
 - 11 human receptors
 - -1 rat (PXR)
 - β-lactamase reporter gene assays except:
 - PXR assays are luciferase reporter gene assays
- p53 Reporter Gene assay
 - β-lactamase reporter gene assay
- Parental cell lines mostly HEK293 (also HeLa and DPX-2)
- 12-15 point concentration-response curves (single replicate)

NCGC: Data Calculations

- United States
 Environmental Protection
 Agency
 - Data normalized to reference compound effect
 - Curves fit to 3- or 4-parameter Hill equation
 - Artifacts removed where obvious fluorescence or cytotoxity detected
 - Required at least 25% efficacy of control compound to calculate AC50
 - AC50 values provided
 - Antagonist format assays challenging due to effects of cytotoxicity
 - LXR assay problematic contaminated with GR reporter line?



■ NCGC00164230-01 ■ NCGC00022570-07



NCGC Assays: What is Being Measured?

- NR assays are ligand-detection assays
- False positives
 - Fluorescent compounds
 - Statistical vs biological significance
 - Gal4_NR-LBD not physiological
 - Cytotoxicity (antagonist format)
 - Operational
- False negatives
 - Fluorescent compounds
 - Cytotoxity
 - Gal4_NR-LBD not physiological
 - Operational
 - Lack of biotransformation



Additional Data Sets To Be Added Soon:

- NHEERL
 - -Zebrafish developmental toxicity (Padilla) Poster
 - Neurite outgrowth and neuronal proliferation (Mundy and Shafer)
 - –ES cell differentiation (Hunter)
- Plasma protein binding and hepatocyte clearance (Thomas)
- PPARα and AhR (NCGC)
- Primary rat hepatocyte HCS Cell Health (Cellumen)