Opportunities and Strategies to Further Reduce Animal Use for Leptospira Vaccine Potency Testing

Angela Walker, DVM, PhD

Center for Veterinary Biologics
Veterinary Services
Animal and Plant Health Inspection Service
United States Department of Agriculture
Ames, Iowa USA







Animal Use



- Regulatory Potency Testing
 - 9CFR 113.101, 113.102, 113.103, 113. 104
- Culture Maintenance









Regulatory Testing

<u>Ideal</u>

- Replacement of hamster vaccinationchallenge model with an in vitro test
 - USDA CVB ELISAs
- CVB Notice 09-16







Regulatory Potency Testing

	9 CFR
Vaccinates	10
Challenge Controls	10
Back-titration	20
Total animals/serial (1 serial tested)	40
Total animals/serial (2 serial tested)	25
Total animals/serial (3 serial tested)	20





<u>Ideal</u>

- * In vitro culture method
 - Virulent
 - Stable
 - Highly prolific
 - "Master Seed" Challenge







Reality

- * Lots of variables affect virulence
 - pH
 - Osmolarity
 - Temperature
 - BSA Source
 - Others



Current Protocols in Microbiology





Strain specificity

- ❖ 250+ pathogenic serovars
 - Variations among leptospiral strains

Photo from The Leptospirosis
Information Center







- In vivo passage
 - Continuous passage of virulent *Leptospira* through hamsters
 - Proposed: Intermittent passage of virulent Leptospira through hamsters
 - Virulence
 - Minimal time to initiate testing
 - Continuous supply of challenge strains





Challenge Maintenance Current Proposed










Culture Maintenance



Parameters

- Quality control of sample
- Speed of Freezing
- Cryopreservative
- **❖** Thaw → Inoculation


Evaluation

- Virulence after short-term and long-term storage
 - Three passages through hamsters
 - Challenge in 10 hamsters





Preliminary Study Example








Culture Maintenance:

Preliminary Results

- Quality control of sample
 - Spirochete Counts
- Speed of Freezing
 - Short-freezes ≥ Gradual Freeze
- Cryopreservative
 - Glycerol > DMSO
 - < 5% Glycerol
- **❖** Thaw → Inoculation



Fluorescent antibody detection of *L. interrogans* serovar *pomona*. Photograph courtesy of Richard Hornsby





Summary

- Regulatory Potency Testing
 - Ideal: ELISAs developed by CVB
 - Reduced Animal Use
 - Shared Controls
 - Re-evaluating back-titrations

Culture Maintenance

- ➤ Ideal: *In vitro* culture
 - Not feasible at this time
- Reduced animal use: Intermittent passage





Thank you

- Geetha Srinivas
- Renee Olsen
- Dave Alt
- Mark Wilson







Questions?





