

Subject: [NTP Web] [NICEATM] INFORMATION OECD 431: 2013

Date: Tuesday, November 26, 2013 9:48:28 AM Eastern Standard Time

From: [REDACTED]

To: NIEHS NICEATM

Name: Miss Francesca Cassata

Organization: Chelab Srl

Email: [REDACTED]

Phone:

Comments: Dear all,
in our laboratory, we are using OECD 431:2013 to convalide RHE inserts.
Skinethic method provide to incubate RHE at 37 °C for 1 hour, while OECD 431 (Exposure time and temperature, pg 27) provide to incubate RHE at room temperature for 1 hour.
What do you suggest? Which is the temperature that you use?

Best regards

Francesca Cassata

Chelab merieux nutriscience