

NTP Approaches to Assessment of Dermal Hypersensitivity

Dori R. Germolec, Ph.D.

Toxicology Branch

NIEHS, DNTP

Board of Scientific Counselors Meeting

December 7-8, 2017





Skin Sensitization

“Allergic Contact Dermatitis”





**NOT TESTED
ON ANIMALS**

Accounts for 10-15% of all occupational disease (*Anderson et al. 2010*)

Major testing requirement for cosmetics, pesticides, industrial chemicals, etc.



Skin Sensitization Process



*Illustration by D. Sailstad



***In Vivo* Tests for Assessment of Dermal Sensitization**



Guinea Pig Maximization Test

- Intradermal and topical sensitization
- Topical challenge
- Measure erythema response 24 - 48 hours post challenge

Buehler

- Topical sensitization with closed patch
- Topical challenge distal to sensitization with closed patch
- Measure erythema response following removal of patch

















Local Lymph Node Assay

- Topical treatment on dorsal surface of the ear
- Inject with radiolabel or fluorochrome
- Measure cell proliferation in the lymph nodes associated with the site of application



U.S. Regulatory Requirements/Considerations

		Reference Animal Method	Classification Criteria
	Pesticides, Industrial chemicals	 LLNA	NS S   Hazard
	Household Products	 LLNA	NS S SS    Potency
	Dermatological Products	 GPMT	   Potency*

Daniel et al. 2017 in preparation

*preference



Accuracy of Animal Tests Against Human Data

GPMT / Buehler



Hazard

~72%

Potency (GHS)

~60%

LLNA



Hazard

72%-82%


Potency (GHS)

54% - 60%

Interagency Coordinating Committee on the Validation of Alternative Methods (ICCVAM). 1999. NIH Publication No. 99-4494
ICCVAM. 2010. NIH Publication No. 11-7709
Urbisch et al. 2015. Reg Tox Pharm 71:337-351.
Hoffmann et al. 2017 in preparation



Reproducibility of LLNA Data



How concordant are multiple LLNA outcomes for a single chemical?

- ~78% for hazard
- ~62% for potency classification






Comparison of LLNA and Human Data

Accuracy:


- **75%** for Hazard (NS/S)
- **60%** for Potency 3-class (NS, Weak/ Moderate, Strong/ Extreme)
- **47%** for Potency 5-class (NS, Weak, Moderate, Strong, Extreme)



Provides a benchmark for comparison with new approaches




Key Events in the Skin Sensitization Process



*Illustration by D. Sailstad




OECD AOP for Skin Sensitization






OECD AOP for Skin Sensitization





Global Skin Sensitization Project

- Objective: analysis of available non-animal approaches
 - OECD submitted case studies
- Collaboration with Cosmetics Europe
 - 128 substance dataset
 - LLNA and human data
 - Curation/generation of *in vitro* data
 - DPRA, KeratinoSens, hCLAT, U-SENS
 - PPRA, SENS-IS (underway)
- Analyze five OECD-submitted defined approaches (i.e., code packages); open source and transparent (R, Python)
- Evaluate performance against the LLNA and human hazard/potency categories



Spectrum of 128 substances
(largely cosmetic ingredients)




In Vitro Models for Assessing Dermal Sensitization

- ***In vitro* assays**
 - **Direct peptide reactivity assay**
 - Assesses the ability of a substance to form a hapten-protein complex
 - **KeratiNoSens**
 - Assesses the ability of a substance to activate cytokines and induce cytoprotective genes in keratinocytes
 - **h-CLAT**
 - Assesses the ability of a substance to activate and mobilize dendritic cells in the skin



DPRA

- Assesses protein reactivity of a test substance
- Uses two heptapeptides
 - One with cysteine (Cys) and one with lysine (Lys) as the reactive center
 - Incubate with test substance and measure disappearance of peptides with HPLC
 - Average depletion (Ave.Lys.Cys) > 6.38% = sensitizer
 - OECD Test Guideline 442C (2015)



Measurements

%Cys depletion


%Lys depletion

%Ave.Lys.Cys depletion

Graphic from EC JRC. 2012. Direct Peptide Reactivity (DPRA) Validation Study Report.



- Assesses the activation of the AKR1C2-ARE element, an indication of keratinocyte activation, in KeratinoSens cells (derived from HaCaT keratinocytes)
 - Caused by electrophilic agents, which tend to be skin sensitizers
 - Measures fold-induction of luciferase activity; induction >1.5-fold in 2/3 experiments = sensitizer
 - OECD Test Guideline 442D (2015)



Measurements


Fold induction

EC1.5

% viability



- Measures 2 cell surface markers, CD86 and CD54, on dendritic cell surrogates (THP-1 cells)
 - Assesses the maturation process of dendritic cells as they transform from antigen processing cells to antigen presenting cells
 - CD86 relative fluorescence intensity (RFI) $\geq 150\%$ and/or CD54 RFI $\geq 200\%$ at any dose, in at least 2/3 experiments, then substance is a sensitizer
 - OECD Test Guideline 442E (2016)





Individual Assays Compared to the LLNA

	hCLAT vs LLNA	DPRA vs LLNA	Keratino vs LLNA																											
	<div>↓</div> <table><tr><td></td><td>NEG</td><td>POS</td></tr><tr><td>NEG</td><td>16</td><td>15</td></tr><tr><td>POS</td><td>17</td><td>79</td></tr></table>		NEG	POS	NEG	16	15	POS	17	79	<div>↓</div> <table><tr><td></td><td>NEG</td><td>POS</td></tr><tr><td>NEG</td><td>22</td><td>32</td></tr><tr><td>POS</td><td>11</td><td>62</td></tr></table>		NEG	POS	NEG	22	32	POS	11	62	<div>↓</div> <table><tr><td></td><td>NEG</td><td>POS</td></tr><tr><td>NEG</td><td>22</td><td>31</td></tr><tr><td>POS</td><td>11</td><td>64</td></tr></table>		NEG	POS	NEG	22	31	POS	11	64
	NEG	POS																												
NEG	16	15																												
POS	17	79																												
	NEG	POS																												
NEG	22	32																												
POS	11	62																												
	NEG	POS																												
NEG	22	31																												
POS	11	64																												
Sensitivity %:	84.0	66.0	67.4																											
Specificity %:	48.5	66.7	66.7																											
Accuracy %:	74.8	66.1	67.2																											
	n=127	n=127	n=128																											



Individual Assays Compared to Human

hCLAT vs Human			DPRA vs Human			Keratino vs Human		
↓			↓			↓		
	NEG	POS		NEG	POS		NEG	POS
NEG	20	11	NEG	29	25	NEG	31	22
POS	19	77	POS	10	63	POS	8	67
Sensitivity %: 87.5			71.6			75.3		
Specificity %: 51.3			74.4			79.5		
Accuracy %: 76.4			72.4			76.6		
≥ LLNA						≥ LLNA		
n=127			n=127			n=128		



Defined Approach Evaluation

- Most non-animal testing strategies evaluated so far perform **better** than the LLNA at predicting human skin sensitization **hazard and potency**
- Combining multiple *in vitro* assays and *in silico* methods or physico chemical properties increases the ability to predict sensitizers



Combining *in vitro* assays and other approaches increases the ability to predict sensitizers

Table 6. LOOCV results for seven highest performing SVM models

No.	Model (Accuracy ^a)	Sensitivity (%)	Specificity (%)	Accuracy (%)
1	DPRA + KeratinoSens + h-CLAT + Toolbox + Lys + Cys + Avg.Lys.Cys + 6 properties (95%)	89	91	89
5	KeratinoSens + h-CLAT + Toolbox + Avg.Lys.Cys + 6 properties (95%)	92	79	88
7	h-CLAT + Toolbox + 6 properties (97%)	85 (84)	94 (48)	88 (75)
8	KeratinoSens + Toolbox + Avg.Lys.Cys + 6 properties (94%)	84 (67)	91 (67)	86 (67)
9	KeratinoSens + h-CLAT + Avg.Lys.Cys + 6 properties (92%)	89	73	84
10	h-CLAT + Toolbox + Avg.Lys.Cys + 6 properties (92%)	90	88	89
11	KeratinoSens + h-CLAT + Toolbox + 6 properties (92%)	89	79	86

Avg.Lys.Cys, average depletion for lysine and cysteine; Cys, average % cysteine; DPRA, direct peptide reactivity assay; h-CLAT, human cell line activation test; LOOCV, leave-one-out cross-validation; Lys, average % lysine depletion; Toolbox, read-across using QSAR Toolbox; SVM, support vector machine.

^aAverage accuracy of the training and test sets for predicting the reference LLNA outcomes.

(Individual assay compared to the LLNA)



Expanding Coverage of Chemical Space

- **Most chemicals used in the validation of non-animal test methods have been cosmetics ingredients**
- NTP is supporting testing of other types of chemicals in three alternative test methods: DPRA, LuSens, hCLAT
 - Expanded chemical space includes: pesticides, agrochemical formulations, dermal excipients, personal care product ingredients, “challenge” chemicals
- Have compiled chemical nominations from multiple ICCVAM agencies
 - EPA: Office of Pesticides, Office of Pollution Prevention and Toxics, Office of Research and Development
 - Consumer Product Safety Commission
 - Food and Drug Administration
 - NTP



Expanding Coverage of Chemical Space

- Total of 266 chemicals nominated
- NTP has procured 135 chemicals for initial testing phase (mostly nominations from the EPA)
- Testing began in late 2017
- Additional testing (~100 chemicals) to follow in mid-2018
- Coordinating with Dow to test formulations already assessed in DPRA and KeratinoSens™ in the hCLAT assay



Expanding Coverage of Chemical Space

- Combine with *in silico* data and physico chemical properties when available
- Evaluate the dataset using methods previously developed by NICEATM (Strickland et al 2016)
- Evaluate predictive performance of non-animal defined approaches submitted to OECD (Kleinstreuer et al. 2018) in comparison to LLNA data
- Characterize applicability domain of in vitro test methods and non-animal defined approaches
- Work with ICCVAM agencies to adopt non-animal defined approaches where appropriate



Acknowledgments

- Nicole Kleinstreuer, Warren Casey (NICEATM)
- Victor Johnson, Michelle Miller (Burleson Research Technologies)
- Bradley Collins (NTP Program Operations)
- David Allen, Judy Strickland, Dan Zang, Mike Paris, Eileen Phillips (ILS)
- Evisabel Craig, Anna Lowit (EPA/OPP)
- Joanna Matheson (CPSC)